

Radius® Spinal System

Surgical Technique

Degenerative



Introduction

The Radius Spinal System incorporates a number of engineering features that makes the system unique. These features include a **one-step (two-stage) locking mechanism**, a floating saddle designed to provide circumferential locking of the rod and a **defined locking position** of the rod into the implant seat.

Most spinal systems use conventional threaded set screws that must be tightened to a predetermined torque limit and are known to occasionally cross thread. The Radius system has a **non-threaded locking mechanism** that requires **no torque measuring device**. **The locking mechanism is designed to reduce the potential for false locking and cross-threading**. This award-winning, non-threaded technology is designed to help increase **the speed, ease, and reliability** of connecting rods to screws.⁺

Acknowledgements

Stryker Spine would like to extend their thanks to the following surgeons for their participation in the development of the Radius System:

- Bruce V. Darden, MD
- Eric B. Laxer, MD
- R. Alden Milam, MD
- Daniel B. Murrey, MD
- Alfred L. Rhyne, MD

Important: The Radius implants and instruments are designed and tested for use only with the Radius Spinal System. This Surgical Technique sets forth detailed, recommended procedures for using the Radius Spinal System implants and instruments. It offers guidance that you should heed but, as with any such technical guide, each surgeon must consider the particular needs of each patient and make appropriate adjustments when necessary and as required. Always refer to the package insert, product label and/or instructions before using any Stryker implant or instrument.

Note: This is intended as a guide only. There are multiple techniques for the insertion of pedicle screws and, as with any surgical procedure, a surgeon must be trained before proceeding.

⁺Edmund Scientific Award for Excellence in Design. Featured in the Design News article on “Fastening innovations speed spine surgery”, 2/26/01.

Table of Contents



Degenerative Surgical Technique

Introduction	02
Key Design Features	04
Indications for Use	06
I. Patient Positioning	06
II. Pedicle Preparation	07
III. Screw Insertion	09
IV. Rod Contouring	12
V. Rod and Cap Insertion	13
VI. Rod Linkage	14
VII. Final Tightening	18
VIII. Cross Connectors	19

Implants	22
Instruments	24

Key Design Features

Ease of Use

An easy, low torque twist locks the construct – it's that simple!

Low Profile

The unique dovetail feature of the helical Wedgelock mechanism allows for a low volume, multi-angle screw and is designed to keep the screw head from spreading during locking.

- Low screw head height is patient friendly
- Small screw head diameter may allow for less facet joint impingement and more room for bone graft
- Short rod run-out is well-suited for L5 – S1 fusions and cross connector placement
- Less overall profile enables better visualization of the anatomy


Defined Locking Position

At the fully locked position, the locking cap is designed to consistently and securely seat the rod to the implant in the optimal position. When in the fully locked position, the laser marked lines on the locking cap are within the laser marked zone on the screw head, giving a visual indicator that the system is fully locked.




Comprehensive top loading spinal rod system consisting of titanium and Vitallium implants and instruments designed for stabilization of the spine.


Large surface area of locking cap saddle designed to **evenly distribute load onto the rod** to help increase sliding strength and reduce risk of notching.




Two-stage (provisionally and fully locked) **locking cap** designed to allow for easy compression and distraction adjustments.



Initial




Provisional




Fully Locked

Locking cap helical wedge ramp facilitates **optimized seating of the rod without** the need for a **torque wrench**.



Telescoping cross connector with angled locking plugs is oriented to allow cross-wound tightening. A multi-angle cross connector is available with a polyaxial head on one side designed to eliminate the need to bend the cross connector to accommodate divergent rods.



The dual lead cortical thread of the **Radius Rapid Screws** is designed to help increase the speed of screw insertion.



Indications for Use

The Radius Spinal System is intended for use in the noncervical spine. When used as an anterior/anterolateral and posterior, noncervical pedicle and non-pedicle fixation system, the Radius Spinal System is intended to provide additional support during fusion using autograft or allograft in skeletally mature patients in the treatment


of the following acute and chronic instabilities or deformities:

- Degenerative disc disease (DDD) (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies);
- Spondylolisthesis;
- Trauma (i.e., fracture or dislocation);
- Spinal stenosis;

- Curvatures (i.e., scoliosis, kyphosis, and/or lordosis);
- Tumor;
- Pseudoarthrosis; and
- Failed previous fusion.


The Radius Spinal System can also be linked to the Xia Titanium Spinal System via the Ø5.5mm to Ø6.0mm Radius rod-to-rod connector.

I. Patient Positioning



Diagnosis is based upon patient history, physical findings, and preoperative radiographic assessment.

The patient can be positioned on the spine table in the prone position. Care should be taken to pad all bony prominences. To facilitate venous drainage, do not compress the abdomen.



Surgical levels may be verified clinically or radiographically. To ensure adequate exposure, the incision is made to extend just beyond the length of the intended fusion.

Use presurgical planning to define the most appropriate implants as well as the optimal location of the implant insertion.

II. Pedicle Preparation


Pedicle Preparation

Remove the small cortical crest with a rongeur or power burr to expose the underlying cancellous bone.

The entry point may be prepared with a **Pedicle Awl**.

Note: It is recommended that you use the **OASYS Awl** (48560010) and **Probe** (48561115) when preparing a pathway for the Radius 4.0mm diameter screws. DO NOT use the **Radius Awl** and **Probe** for this screw diameter.

Pedicle Awls





A pathway is then opened with a **Probe**. The **Probe** should contact the bone at all times. The correct rotational insertion of the instrument will allow the **Probe** to follow a path of least resistance without violating the pedicle walls.

In the event that resistance is felt, the entry point and trajectory should be reevaluated.



The **Probe** is laser etched with 10mm intervals to indicate the depth in which the **Probe** has been inserted as well as to help determine proper screw length.

Probes




The **Probes** have been colored gold up to 40mm to allow easy visualization of 40mm depth, which represents the most common screw length.

Pedicle Testers



Check the prepared pathway with a **Pedicle Tester** to verify that all walls of the pedicle have not been violated and that cancellous bone is felt at the distal end of the path.

II. Pedicle Preparation



The appropriate **Tap** may be used to prepare the pedicle canal when the surgeon is having difficulty starting the screw. The **Tap** sizes for degenerative surgery are Ø3.5mm, Ø4.25mm, Ø5.25mm, Ø6.25mm, and Ø7.25mm. Radius also has iliac fixation with tap sizes of Ø8.25mm and Ø9.25mm.

Note: Use the **Taps** when preparing the pedicle canal for Standard and Multi-Angle Screws. Use the **Rapid Taps** when using Radius Rapid Multi-Angle Screws. The **Rapid Taps** for degenerative surgery come in diameters of 5.25mm, 6.25mm, and 7.25mm. Rapid Taps for iliac fixation come in diameters of 8.25mm and 9.25mm.

The **Taps** are laser etched with 10mm intervals to indicate the depth to which they have been inserted as well as to help determine proper screw length. The **Taps** are 0.5mm smaller than the corresponding screw and feature a spiral cutting flute to minimize toggle during tapping.

The **Taps** have been colored gold to 40mm to allow easy visualization of 40mm depth, which represents the most common screw length.

The **Drill Sleeve** serves as both a soft tissue protector and a depth indicator and is intended to be used with the **Taps**.



Note: The **Drill Sleeve** color rings denote which **Tap** sizes fit through the Sleeve.

Tip: If soft tissue unintentionally redirects the **Tap** laterally toward the pedicle, the **Drill Sleeve** can be used to lever against the soft tissue without bending the **Tap**.

Tip: It is recommended that you use the **Drill Sleeve** for the Ø3.5mm **Tap**.

III. Screw Insertion



Both the **Standard Screw Inserter** and the **Multi-Angle Screw Inserter** are designed to provide a very rigid connection with the **Standard Screws** and **Multi-Angle Screws**.

With the pedicle pathways prepared and proper screw length and diameter determined, the screw is ready for insertion.

Standard Screws


To assemble the Standard Screw Inserter:

1. Insert the inner shaft through the body of the **Standard Screw Inserter**.
2. Insert the ratchet down the shaft of the **Standard Screw Inserter**. Verify that the ratchet is bottomed out.
3. Connect to the desired handle.
4. Ensure that the **Standard Screw Inserter** is fully unlocked.
5. Align the **Standard Screw Inserter** shaft with the internal geometry of the standard screw head while holding the bone screw.
6. Rotate the dial on the **Standard Screw Inserter** clockwise to firmly seat the screw onto the **Standard Screw Inserter**.



To insert Standard Screw:

1. Ensure that the **Standard Screw Inserter** is fully unlocked.
2. Drop the **Standard Screw Inserter** into the head of the screw.
3. Rotate the dial on the **Standard Screw Inserter** clockwise to firmly lock the screw onto the **Standard Screw Inserter**.
4. Once the screw is inserted into the pedicle, push the gold button and rotate the dial on the **Standard Screw Inserter** counter-clockwise until the **Standard Screw Inserter** is loose from the screw.
5. Remove the **Standard Screw Inserter** from the screw.



Standard Screw Inserter

Tip: The **Standard Screw Inserter** can also be used to hold hooks for insertion.

Note: Use the **Rocker** for removal of the **Standard Screws**.

III. Screw Insertion



Multi-Angle Screws

To assemble the Multi-Angle Screw Inserter:

1. Insert the inner shaft through the body of the Multi-Angle Screw Inserter.
2. Insert the ratchet down the shaft of the Multi-Angle Screw Inserter. Verify that the ratchet is bottomed out.
3. Connect to the desired handle.
4. Ensure that the Multi-Angle Screw Inserter is fully unlocked.
5. Align the tabs on the Multi-Angle Screw Inserter shaft with the external quad on the screw head while holding the bone screw.
6. Rotate the dial on the Multi-Angle Screw Inserter clockwise to firmly seat the screw onto the Multi-Angle Screw Inserter.


To insert Multi-Angle Screw:

For Multi-Angle Screws use the Multi-Angle Screw Inserter.

To use the Multi-Angle Screw Inserter:

1. Ensure that the Multi-Angle Screw Inserter is fully unlocked.
2. Align the tabs on the Multi-Angle Screw Inserter shaft with the external quad on the screw head while holding the bone screw.
3. Rotate the dial on the Multi-Angle Screw Inserter clockwise to firmly seat the screw onto the Multi-Angle Screw Inserter.


Note: The Multi-Angle Screw Inserter is used for both the Multi-Angle Screws and the Rapid Multi-Angle Screws.



Multi-Angle Screw Inserter

4. Once the screw is completely inserted into the pedicle, depress the gold button and rotate the dial on the **Multi-Angle Screw Inserter** counter-clockwise until the **Multi-Angle Screw Inserter** is loose from the screw.
5. Remove the **Multi-Angle Screw Inserter** from the head of the screw.

Tip: When using either the **Standard Screw Inserter** or the **Multi-Angle Screw Inserter** there are 4 choices of handles which can be used: **Cannulated T-Handle**, **Cannulated Ratcheting T-Handle**, **Cannulated Standard Handle** or **Cannulated Standard Ratcheting Handle**.




Height Adjustment

The **Multi-Angle Screw** can more easily be adjusted once you remove the **Screw Inserter**, using the **Multi-Angle Adjustment Tool**. Simply attach the **Adjustment Tool** to the inner cruciform and turn counter-clockwise to raise the screw. The **Adjustment Tool** is able to engage the **Multi-Angle Screw** at $\pm 30^\circ$ in any direction.



IV. Rod Contouring



Once all screws are inserted, use the appropriate precut and/or prebent rods.

Rod bending may be needed to fit the desired spinal contours.

Note: Care should be taken to not make extreme bends as that can cause stress, concentration, and notching of the rod.



Bending can be performed with the **Rod Bender**. To contour the rod, a series of small incremental adjustments will bend the rod gradually and provide even stress distribution along the rod.



Rod bending can also be performed using the bending iron feature on the **In-Situ Benders**. Insert the rod into the holes on each bender and gently pull upwards on one bender to contour the rod.

Tip: The **Rod Template** (03710620) can be used to measure the necessary curvature of the rod prior to contouring.

V. Rod and Cap Insertion

Rod Insertion

Once the rod is bent to the desired contour, a **Rod Insertion Forceps** or **Rod Insertion Scissor** can be used to place the rod into the grooves of the implant.



Cap Insertion



1. Assemble the **Locking Cap** onto the **Initial Insert** by using thumb pressure or pushing firmly onto a **Locking Cap** in the locking cap caddy (recommended).

Note: The **Initial Insert** Tube can be used to guide the **Locking Cap** into the screw head.

Tip: When using the **Initial Insert**, make sure the laser lines on the **Initial Insert** line up with the laser lines on the **Locking Cap**.

2. Place the **Locking Cap** into the screw head, apply downward pressure to ensure the **Locking Cap** is firmly seated against the rod, then rotate clockwise until the cap snaps into the provisional locked position. This position is indicated when the laser markings on the cap and the **Initial Insert** are parallel with the rod. The **Locking Cap** is now provisionally seated, the rod is blocked and the **Initial Insert** can be removed.

Tip: Do not turn the **Locking Cap** past the provisional position with the **Initial Insert**. This must be done with the **Final Driver** (see the Final Tightening section).



VI. Rod Linkage

Rod Pusher



The Radius System offers three options for linking the rod to the spine:


Option 1:

Rod Pusher and Initial Inserter

The **Rod Pusher** is designed to push the rod down in-situ and is helpful for small adjustments.

While utilizing the **Rod Pusher**, perform the cap insertion technique.

Note: Do not perform final tightening of the **Locking Cap** with the **Initial Inserter**. It will not be possible to completely lock the **Locking Cap**.




Option 2:

Rocker and Initial Inserter

The **Rocker** is used when the rod is slightly proud with respect to the seat of the implant.

The **Rocker** is designed to easily slide into the lateral holes on the implant head. Lever the **Rocker** back. This forces the rod into the head of the implant. Insert the **Locking Cap** with the **Initial Inserter** when the rod is fully seated in the head of the implant.



Rocker




VI. Rod Linkage


Option 3:

Persuader and Capspin

Use the **Rod Persuader** when additional force is needed to bring the rod to the implant. The **Rod Persuader** is designed to use the **Locking Cap** to persuade the rod directly into the screw head and automatically load the **Locking Cap** into the provisional locked position.




1. Ensure that the lower Bar is in the neutral position. Load the **Locking Cap** onto the **Capspin Assembly** of the **Rod Persuader** by pushing firmly on the **Locking Cap** in the **Locking Cap Magazine** (recommended).



Tip: Rotate the **Capspin Assembly** counter-clockwise, using the **Upper T-Handle** to avoid allowing the **Locking Cap** to contact the rod early, which would prevent the **Rod Persuader** from attaching to the screw head.

2. With the **Locking Cap** assembled, place the **Rod Persuader** onto the screw head by clamping the forceps into the lateral holes on the implant head.



VI. Rod Linkage



3. Turn the Upper T-Handle clockwise until the rod is fully seated. The bottom T-Bar will rotate and click into the provisionally locked position indicating that the rod is completely captured and the locking cap is provisionally locked.

Tip: If the T-Bar does not turn, hand tighten or loosen and try again.



4. Once the Bar twists to provisional lock, turn the T-Handle 3 or 4 rotations counter-clockwise. Then lift the Rack and remove the Rod Persuader.

Tip: Make sure the T-Bar does not rotate backwards.

Tip: The rod cannot be linked to the screws or the hooks if the rod has a sharp, acute bend at the point of linkage.


Tip: The Rod Persuader is not designed to bend the rod.

VI. Rod Linkage

Note: In the event that the rod is forced down while tightening be sure that the **Locking Cap** is fully engaged in the screw head.

Extra caution is advised when:

1. The rod is not horizontally placed into the screw head.
2. The rod is high in the screw head.
3. An acute convex or concave bend is contoured into the rod.






Compression and Distraction

To compress or distract:


1. Final tighten one of the screws.
2. Place the **Counter Torque Tube** instrument over the provisionally locked adjacent screw.
3. Compress or distract using the **Compressor** or **Distractor**.

Tip: Compression and distraction using the under rod instruments depends on how proud the screw is since the instruments may not fit under the head of the screw.



4. After compression or distraction has occurred, place the **Final Driver** through the **Counter Torque Tube** to tighten.

VII. Final Tightening



Once the correction procedures have been carried out and the spine is fixed in a satisfactory position, the final tightening of the **Locking Caps** is done by utilizing the **Final Driver** and the **Counter Torque Tube**.

1. Place the **Counter Torque Tube** around the screw head and over the rod.
2. Place the **Final Driver** into the star on the **Locking Cap** by first aligning the lines on the **Final Driver** shaft with the lines on the **Locking Cap**.
3. Twist the **Final Driver** until the handle of the **Final Driver** is perpendicular to the handle of the **Counter Torque Tube** or until the **Final Driver** no longer turns. The laser marked lines on the **Locking Cap** should be within the laser marked zone on the screw head.

Tip: It is important to position the **Counter Torque Tube** before the **Final Driver** is inserted to minimize the risk of instrument slippage.

Tip: For final locking, the line on the **Final Driver** should be centered in the distal window.

Note: The **Counter Torque Tube** must be used for final tightening. The **Counter Torque Tube** is designed to perform three important functions:

1. To align the **Final Driver** with the tightening axis.
2. To maximize the torque needed to lock the implant assembly.
3. To hold the construct in place during final tightening.

Note: For removal, reverse the final tightening sequence.

VIII. Cross Connectors

Cross Connectors

Cross Connectors are recommended for increased rotational stability of the construct.

Use the Cross Connector Caliper to determine the appropriate size cross connector.

Choose the appropriate Variable or Fixed Cross Connector length and place it onto the rod.

Note: The Variable Multi-Angle Curved Cross Connector is designed to be used when only a portion of the spinous process is removed.

Fixed Cross Connector
(22mm - 36mm)

Variable Cross Connector
(Short, Med, Long; 35mm-86mm)

Variable Multi-Angle Cross Connector
(Med, Long; 46mm-75mm)

Variable Multi-Angle Curved Cross Connector
(XS, S, Med, Long, XL, XXL; 35mm-104mm)

Accommodating Divergent Rods

Two options are available to accommodate divergent rods. The Multi-Angle Cross Connector features a variable head on one side which is designed to allow the connector to be placed on divergent rods without bending.

If using the Variable Cross Connector or the Fixed Cross Connector, bending may be performed to accommodate divergent rods.

The set of Cross Connector Benders has been etched to indicate which slots will accommodate Fixed or Variable Cross Connectors.

Bending Fixed Cross Connectors

1. To bend a Fixed Cross Connector, place either Fixed Cross Connector Clamp into the Fixed Cross Connector Slot of the Cross Connector Bender.
2. Place the other Clamp of the Fixed Cross Connector into the Fixed Nest Slot built into the top of the Cross Connector J-Tip Bender.

Tip: The Fixed Cross Connector may be bent up to 22°. Overbending will result in immediate implant fatigue and possible breakage.




VIII. Cross Connectors



Cross Connector Bender Left



Cross Connector Bender Right



Bending Variable Cross Connectors

1. To bend the **Variable Cross Connector**, place the **Variable Cross Connector Clamp** into the “Variable” Slot of the **Cross Connector Bender**.

Tip: To alleviate confusion, the Slot on the **Cross Connector Bender** will only accommodate the Clamp on the rod side of the **Variable Cross Connector**.

2. Grasp the internal bar of the **Variable Cross Connector** with the “Variable” Nest Slot built into the top of the **Cross Connector Bender**. Bend the internal bar.

Tip: It is recommended that the bend be introduced into the rod as close to the Clamp as possible so as not to interfere with the locking procedure.



Cross Connector Bender Left




Cross Connector Bender Right



VIII. Cross Connectors

Cross Connector Tightening

1. Use the **Cross Connector Plug Driver** to turn the Locking Plug in either direction until the line on the Locking Plug is aligned with the laser marked line on the Clamp. Using the cross connector holder, place the connector on the rods and press downward until the clamp engages the rod.



2. Stabilize the **Cross Connector** using a finger in the middle of the **Cross Connector** and tighten the Locking Plug on the opposite side.




3. If using a **Variable Cross Connector**, use the **Cross Connector Closer** to squeeze the Locking Collet to completely lock the **Variable Cross Connector**.

Tip: The gold jaw of the **Cross Connector Closer** is placed against the gold collet of the **Cross Connector** during tightening.

Tip: It is easier to lock the Fixed Head of the **Variable Cross Connector** first.

Apply bone graft to the fusion site and close in the usual manner.



Implants


Reference Number

Description



486610000

Locking Cap



4866113 (15) to (55) Ø4.00mm

4866114 (25) to (55) Ø4.75mm

4866115 (25) to (65) Ø5.75mm

4866116 (25) to (65) Ø6.75mm

4866117 (25) to (65) Ø7.75mm

Multi-Angle Screws



4866185 (25) to (65) Ø5.75mm

4866186 (25) to (00) Ø6.75mm

4866187 (25) to (00) Ø7.75mm

Multi-Angle, Rapid Screws



4866103 (15) to (55) Ø4.00mm


4866104 (25) to (55) Ø4.75mm

4866105 (25) to (65) Ø5.75mm

4866106 (25) to (65) Ø6.75mm


4866107 (25) to (65) Ø7.75mm

Screws, Standard



486614100 - Short
486614110 - Medium
486614120 - Long

Variable Cross Connectors
(35mm-86mm)



4866140(22) to (36)

Cross Connectors
(22mm-36mm)



486614170 - Medium
486614180 - Long

Variable Multi-Angle
Cross Connectors
(46mm-75mm)



486614128 - Extra Small
486614129 - Small
486614130 - Medium
486614140 - Long
486614150 - Extra Long
486614160 - Extra Extra Long

Variable Multi-Angle
Curved Cross Connectors
(35mm-104mm)

Reference Number

Description

486615030	Ø5.5 x 30mm Titanium Rad Rod
486615035	Ø5.5 x 35mm Titanium Rad Rod
486615040	Ø5.5 x 40mm Titanium Rad Rod
486615045	Ø5.5 x 45mm Titanium Rad Rod
486615050	Ø5.5 x 50mm Titanium Rad Rod
486615060	Ø5.5 x 60mm Titanium Rad Rod
486615070	Ø5.5 x 70mm Titanium Rad Rod
486615080	Ø5.5 x 80mm Titanium Rad Rod
486615090	Ø5.5 x 90mm Titanium Rad Rod
486615000	Ø5.5 x 100mm Titanium Rad Rod
486615010	Ø5.5 x 110mm Titanium Rad Rod
486615020	Ø5.5 x 120mm Titanium Rad Rod



486615550	Ø5.5 x 50mm Titanium Max Rad Rod
486615560	Ø5.5 x 60mm Titanium Max Rad Rod
486615570	Ø5.5 x 70mm Titanium Max Rad Rod
486615580	Ø5.5 x 80mm Titanium Max Rad Rod
486615590	Ø5.5 x 90mm Titanium Max Rad Rod
486615500	Ø5.5 x 100mm Titanium Max Rad Rod
486615510	Ø5.5 x 110mm Titanium Max Rad Rod
486615520	Ø5.5 x 120mm Titanium Max Rad Rod



486613003	Ø5.5 x 30mm Titanium Spinal Rod, without Hex
486613053	Ø5.5 x 35mm Titanium Spinal Rod, without Hex
486613004	Ø5.5 x 40mm Titanium Spinal Rod, without Hex
486613054	Ø5.5 x 45mm Titanium Spinal Rod, without Hex
486613005	Ø5.5 x 50mm Titanium Spinal Rod, without Hex
486613006	Ø5.5 x 60mm Titanium Spinal Rod, without Hex
486613007	Ø5.5 x 70mm Titanium Spinal Rod, without Hex
486613008	Ø5.5 x 80mm Titanium Spinal Rod, without Hex
486613009	Ø5.5 x 90mm Titanium Spinal Rod, without Hex
486613000	Ø5.5 x 100mm Titanium Spinal Rod, without Hex
486613010	Ø5.5 x 110mm Titanium Spinal Rod, without Hex
486613020	Ø5.5 x 120mm Titanium Spinal Rod, without Hex



486613602	Ø5.5mm x 600mm Vitallium Spinal Rod, with Hex
486613601	Ø5.5mm x 600mm Vitallium Spinal Rod, without Hex
486613242	Ø5.5mm x 240mm Vitallium Spinal Rod, with Hex
486613241	Ø5.5mm x 240mm Vitallium Spinal Rod, without Hex

















48661151 (60) to (20)	Ø5.5mm Vitallium Rad Rod, with Hex (60mm-120mm)
-----------------------	--














48661152 (60) to (20)	Ø5.5mm Titanium Rad Rod, with Hex (60mm-120mm)
-----------------------	---















Instruments

	Reference Number	Description
	486619400	Screw Insertion Tray
	486619000	Pedicule Awl
	486619001	Pedicule Awl Conical
	486619002	Pedicule Awl Conical with Stop
	486619030	Sharp Probe Straight
	486619035	Sharp Probe Curved
	486619005	Blunt Probe Straight
	486619010	Blunt Probe Curved
	486619020	Blunt Probe Shaft
	486619025	Blunt Probe Shaft Grooved
	486619085	Probe Shaft T-Handle
	486619041	Pedicule Tester Malleable
	486619043	Pedicule Tester Flexible
	486619040	Pedicule Tester

Reference Number	Description	
486619050	Drill Sleeve	<div> <div> Ø4.00mm Ø4.75mm Ø5.75mm Ø6.75mm Ø7.75mm </div>  </div>
486619052 486619055 486619056 486619057 486619058	Ø3.50mm Tap Ø4.25mm Tap Ø5.25mm Tap Ø6.25mm Tap Ø7.25mm Tap	
486619060 486619061 486619062	Ø5.25mm Rapid Tap Ø6.25mm Rapid Tap Ø7.25mm Rapid Tap	
486619095	Multi-Angle Screw Inserter	
486619110	Standard Screw Inserter	
486619100	Multi-Angle Adjustment Tool	
486619080	Cannulated T-Handle	
486619070	Cannulated Ratcheting T-Handle	
486619075	Cannulated Standard Handle	
486619065	Cannulated Standard Ratcheting Handle	
486619120	3.5mm Hex Shaft	
486619410	Fixation Tray	
486619480	Rapid Fixation Tray	

Instruments

Instruments	Reference Number	Description
	486619170	Rod Bender
	486619145	Rod Insertion Scissor
	486619146	Rod Insertion Forceps
	486619180	Rod Pusher
	486619185	Rocker
	486619186	Locking Cap Magazine
	486619190	Rod Persuader
	486619210	Capspin Assembly
	486619215	Initial Inserter
	486619220	Initial Inserter Tube
	486619380	Under Rod Distractor
	486619390	Under Rod Compressor

Reference Number	Description	
48026100	Compressor	
48026000	Distractor	
486619395	Scoliosis Compressor	
486619385	Scoliosis Distractor	
486619235	Counter Torque Tube	
486619240	Final Driver	
486619431	Cross Connector Tray	
486619280	Cross Connector Caliper	
48040247	Cross Connector Holder	
486619255	Cross Connector Bender Left	
486619260	Cross Connector Bender Right	
486619250	Cross Connector Plug Driver	
486619245	Cross Connector Closer	

Indications

The Radius Spinal System is intended for use in the noncervical spine. When used as an anterior/anterolateral and posterior, noncervical pedicle and non-pedicle fixation system, the Radius Spinal System is intended to provide additional support during fusion using autograft or allograft in skeletally mature patients in the treatment of the following acute and chronic instabilities or deformities:

- Degenerative disc disease (DDD) (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies);
- Spondylolisthesis;
- Trauma (i.e., fracture or dislocation);
- Spinal stenosis;
- Curvature (i.e., scoliosis, kyphosis, and/or lordosis);
- Tumor;
- Pseudoarthrosis; and
- Failed previous fusion.

The Radius Spinal System can also be linked to the Xia Titanium Spinal System via the Ø5.5mm to Ø6.0mm Radius rod-to-rod connector.

Contraindications

Contraindications may be relative or absolute. The choice of a particular device must be carefully weighed against the patient's overall evaluation. Circumstances listed below may reduce the chances of a successful outcome:

- Any abnormality present which affects the normal process of bone remodeling including, but not limited to, severe osteoporosis involving the spine, bone absorption, osteopenia, primary or metastatic tumors involving the spine, active infection at the site or certain metabolic disorders affecting osteogenesis.
- Insufficient quality or quantity of bone which would inhibit rigid device fixation.
- Previous history of infection.
- Excessive local inflammation.
- Open wounds.
- Any neuromuscular deficit which places an unusually heavy load on the device during the healing period.
- Obesity. An overweight or obese patient can produce loads on the spinal system which can lead to failure of the fixation of the device or to failure of the device itself.
- Patients having inadequate tissue coverage of the operative site.
- Pregnancy.
- A condition of senility, mental illness, or substance abuse. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the use of the implant, leading to failure, or other complications.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made prior to material selection or implantation.
- Other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of tumors, congenital abnormalities, elevation of sedimentation rate unexplained by other diseases, elevation of white blood cell count (WBC), or marked left shift in the WBC differential count.

These contraindications can be relative or absolute and must be taken into account by the physician when making a decision. The above list is not exhaustive.

Cautions and Warnings

CAUTIONS (U.S.A.)

Federal law (U.S.A) restricts this device to sale by or on the order of a licensed physician.

The implantation of pedicle screw spinal systems must be performed only by experienced spinal surgeons with specific training in the use of this pedicle screw spinal system because this is a technically demanding procedure presenting a risk of serious injury to the patient.

Based on the fatigue testing results, the physician/surgeon must consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc., which may impact on the performance of the system.

WARNING (U.S.A.)

The safety and effectiveness of pedicle screw spinal systems have been established only for spinal conditions with significant mechanical instability or deformity requiring fusion with instrumentation. These conditions are significant mechanical instability or deformity of the thoracic, lumbar, and sacral spine secondary to spondylolisthesis (grades 3 and 4) of the L5-S1 vertebrae, degenerative spondylolisthesis with objective evidence of neurological impairment, fracture, dislocation, scoliosis, kyphosis, spinal tumor, and failed previous fusion (pseudoarthrosis). The safety and effectiveness of these devices for any other conditions are unknown.

Radius implant components have not been tested for heating or migration in MR environment.

Removal of Implants

These implants are temporary internal fixation devices designed to stabilize the operative site during the normal healing process. After healing occurs, these devices serve no functional purpose and can be removed. Removal may also be recommended in other cases, such as:

- Corrosion with a painful reaction
- Migration of the implant, with subsequent pain and/or neurological, articular or soft tissue lesions
- Pain or abnormal sensations due to the presence of the implants
- Infection or inflammatory reactions
- Reduction in bone density due to the different distribution of mechanical and physiological stresses and strains
- Failure or mobilization of the implant

Standard ancillaries provided by Stryker Spine can be used to remove the implants. Any decision by a physician to remove the internal fixation device must take into consideration such factors as the risk to the patient of the additional surgical procedure as well as the difficulty of removal. Removal of an unloosened spinal screw may require the use of special instruments to disrupt the interface at the implant surface. This technique may require practice in the laboratory before being attempted clinically. Implant removal must be followed by adequate postoperative management to avoid fracture or re-fracture. Removal of the implant after fracture healing is highly recommended. Metallic implants can loosen, bend, fracture, corrode, migrate, cause pain or stress shield bone.

Information for Patients

The surgeon must discuss all physical and psychological limitations inherent to the use of the device with the patient. This includes the rehabilitation regimen, physical therapy, and wearing an appropriate orthosis as prescribed by the physician. Particular discussion should be directed to the issues of premature weightbearing, activity levels, and the necessity for periodic medical follow-up.

The patient must be warned of the surgical risks and made aware of possible adverse effects. The patient must be warned that the device cannot and does not replicate the flexibility, strength, reliability or durability of normal healthy bone, that the implant can break or become damaged as a result of strenuous activity or trauma, and that the device may need to be replaced in the future. If the patient is involved in an occupation or activity which applies inordinate stress upon the implant (e.g., substantial walking, running, lifting, or muscle strain) he/she should be warned that resultant forces can cause failure of the device. Patients who smoke have been shown to have an increased incidence of non-unions. Such patients should be advised of this fact and warned of the potential consequences. For diseased patients with degenerative disease, the progression of degenerative disease may be so advanced at the time of implantation that it may substantially decrease the expected useful life of the appliance. In such cases, orthopaedic devices may be considered only as a delaying technique or to provide temporary relief.

REUSE

Never reuse or reimplant spinal surgical implants. These could become contaminated resulting in infection. In addition, even though the device appears undamaged, it may have small defects which could compromise structural integrity reducing its service life and/or leading to patient injury.

It is recommended to verify that the instruments are in good condition and operating order prior to use during surgery.

Notes

Notes

Reconstructive

Hips
Knees
Trauma & Extremities
Foot & Ankle
Joint Preservation
Orthobiologics & Biosurgery

MedSurg


Power Tools & Surgical Accessories
Computer Assisted Surgery
Endoscopic Surgical Solutions
Integrated Communications
Beds, Stretchers & EMS
Reprocessing & Remanufacturing

Neurotechnology & Spine

Craniomaxillofacial
Interventional Spine
Neurosurgical, Spine & ENT
Neurovascular
Spinal Implants



Stryker Spine
2 Pearl Court
Allendale, NJ 07401-1677 USA
t: 201-760-8000
www.stryker.com



Stryker France
ZAC – Avenue de Satolas Green
Pusignan 69330
France

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: OASYS, Radius, Stryker, Vitallium, Wedgelock, Xia. All other trademarks are trademarks of their respective owners or holders.