



stryker®

# Apex & HA Apex Pins

## Pin Fixation System

- Half Pins, Transfixing Pins
- HA Coated Half Pins for long term fixation
- Instruments



---

This publication sets forth detailed recommended procedures for using Stryker Osteosynthesis devices and instruments.

It offers guidance that you should heed, but, as with any such technical guide, each surgeon must consider the particular needs of each patient and make appropriate adjustments when and as required.

A workshop training is recommended prior to first surgery.

All non-sterile devices must be cleaned and sterilized before use. Follow the instructions provided in our Instructions for Cleaning, Sterilization, Inspection and Maintenance (L24002000). Multi-component instruments must be disassembled for cleaning. Please refer to the corresponding assembly/disassembly instructions.

See package insert (V15011, V15013, V15034) for a complete list of potential adverse effects, contraindications, warnings and precautions. The surgeon must discuss all relevant risks, including the finite lifetime of the device, with the patient, when necessary.

Stryker Osteosynthesis systems have not been evaluated for safety in MR environment and have not been tested for heating or migration in the MR environment, unless specified otherwise in the product labeling.

**Warning:**

**Fixation Screws:**

**Stryker Osteosynthesis bone screws are not approved or intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.**

# Contents

	<b>Page</b>
<b>1. Introduction</b>	<b>4</b>
<b>2. Indications, Features and Benefits</b>	<b>5</b>
<b>3. Implants</b>	<b>7</b>
<b>4. Instruments</b>	<b>8</b>
<b>5. Technical Details</b>	
Instructions for Instrument Use	10
<b>6. Technical Details</b>	
Pin Insertion Guidelines	11
HA Apex Coated Pin Insertion Guidelines	13
<hr/>	
<b>7. Ordering Information – Implants</b>	<b>14</b>
<b>8. Ordering Information – Instruments</b>	<b>19</b>

# Introduction

The Apex Pin line has been a success for more than 20 years.

The self-drilling pin technology was introduced in 1987, and Apex Pins continue to be widely used each day throughout the world<sup>1</sup>.

Solid pin fixation and appropriate pin size is essential for effective external fixator frames.<sup>2</sup> A well designed pin can help to improve treatments and may reduce the risk of complications.<sup>3</sup>

The Apex Pin range offers a wide selection of pins in various lengths and diameters. The Stainless Steel and Titanium Self-Drilling/Self-Tapping Pins offer a one-step insertion where pre-drilling is not required. The self-tapping Stainless Steel Blunt Pins require pre-drilling. Stainless Steel Cancellous Pins include a cancellous thread, designed to provide grip in cancellous bone and require pre-drilling as well. Self-Drilling Stainless Steel Transfixing Pins are available threaded or smooth and are indicated for bilateral frame constructs.

Stryker also offers a range of Hydroxylapatite (HA) coated Apex Pins.

Hydroxylapatite has a long clinical history and is chemically similar to the mineral components of bone. HA is one of the few materials that supports bone ingrowth and osteointegration<sup>4</sup>.

The HA coated Apex Pins are Stainless Steel and are available in Self Drilling/Self Tapping and Blunt. They are offered in 5mm and 6mm diameters with numerous thread and shank lengths to meet the different needs of each application.

Due to Stryker's coating technology the Self Drilling/Self Tapping HA Apex Pins provide excellent cutting and drilling capabilities without compromising the HA coating.

This wide selection of options allows you to choose the most appropriate Apex Pin for your needs.



1. Encompass Sales Data

2. Eric Ledet, PhD., Director of the Orthopaedics Research Library, Albany Medical College; Biomechanical Factors in External Fixation and Hybrid External Fixation. Stryker White Paper 2004. Literature# LSA48

3. Wikenheiser MA, Market MD, Lewallen DG, et al. Thermal response and torque resistance of five cortical half pins under simulated insertion technique. J.Orthop Res 1995; 13: 615-619

4. Moroni, Antonio; Orienti, Luca; Stea, Susanna; Visentin, Manuela. Improvement of the Bone-Pin Interface with Hydroxyapatite Coating: An In Vivo Long-Term Experimental Study. Journal of Orthopaedic Trauma. 10(4):236-242, May 1996



# Indications

Apex Pins are used in conjunction with diverse Stryker External Fixation Systems including the Hoffmann product line, Triax or Tenxor. For Indications please refer to the specific Operative Technique of the External Fixation System. For Contraindications and Precautions always refer to the Instructions for Use in the Package Insert.

## Self-Drilling / Self-Tapping Apex Pin

### Pin Design

When using the self-drilling Apex Pin, pre-drilling is not necessary. The Self Drilling/Self Tapping pins' cutting geometry potentially allows for reduced insertion temperature<sup>3</sup>.



A double helical flute creates a homogeneous thread profile that transports bone chips out of the drill hole.

The thread design maximizes contact with the bone and controls stress distribution on the pin/bone interface by optimizing radial tension<sup>3</sup>.

**Self-Tapping section**

**Self-Drilling tip**

3. Wikenheiser MA, Market MD, Lewallen DG, et al. Thermal response and torque resistance of five cortical half pins under simulated insertion technique. J.Orthop Res 1995; 13; 615-619

# Implants

## HA Coated Apex Pin

Hydroxylapatite coating is one of the few materials that supports bone ingrowth and osteointegration<sup>4</sup>. The HA coated Self Drilling/Self Tapping and Blunt Apex Pins come in a range of diameters, lengths and thread lengths to provide you with a broad range of options.

Blunt HA coated Apex Pins require pre-drilling.




Although HA coated, this Pin allows for a one step procedure due to its self-drilling and self-cutting design technology.



HA coating along the full length of the threads provides an excellent pin/bone interface in both cortices.



The double helical flute creates a homogeneous thread profile that transports bone chips out of the drill hole for additional improvement of the pin/bone interface<sup>3</sup>.



Pre-drilling is required for this pin.

The cylindrical pin design helps increase bone purchase and pull out resistance. Since the threads don't taper, the HA Apex Pins can be backed out without compromising fixation.



3. Wikenheiser MA, Market MD, Lewallen DG, et al. Thermal response and torque resistance of five cortical half pins under simulated insertion technique. J. Orthop Res 1995; 13: 615-619

4. Moroni, Antonio; Orienti, Luca; Stea, Susanna; Visentin, Manuela. Improvement of the Bone-Pin Interface with Hydroxyapatite Coating: An In Vivo Long-Term Experimental Study. Journal of Orthopaedic Trauma. 10(4):236-242, May 1996

# Implants

## Apex Pin Range


### **Self-Drilling / Self-Tapping Pin**

Stainless Steel and Titanium Self-Drilling/Self-Tapping Pins allow a one-step procedure due to their self-drilling and cutting design technology.




### **Blunt / Self-Tapping Pin**

Stainless Steel Blunt Pins require pre-drilling.



### **HA Coated Pin**

The range of HA coated Apex Pins offers a wide variety of Stainless Steel Self Drilling/Self Tapping and Blunt Pins in numerous lengths, thread lengths and diameters to meet the different needs of each application.




### **Cancellous Pin**

Stainless Steel Cancellous Pins are designed for a strong grip in cancellous bone. The cancellous thread provides increased contact area between the cancellous bone and the pin. This pin is blunt and requires pre-drilling.



### **Transfixing Pin**

Self-Drilling Transfixing Pins are available threaded or smooth and are indicated for bilateral frame constructs.




# Instruments

## Apex Instruments

### Drill Brace

The Drill Brace is designed for manual pin insertion for better control and reduced insertion temperature.

It provides integrated attachments for 3mm & 4mm and 5mm & 6mm pins. Simply by changing the Drill Handle from one end to the other you gain access to the different attachments.




The spring mechanism is designed to provide a secure pin connection during insertion. The curvature of the brace gives a smooth insertion feel.

### Drill Guide

The Drill Guide is designed for simplified parallel pin insertion. The color coded Drill Guide Blocks provide the correct distance for the various pin clamps of the Stryker's External Fixation Systems (e.g. Hoffmann 3, Hoffmann II MRI, Hoffmann II Compact MRI, and Monotube Triax and Tenxor systems.) The color coding matches the colors of the various systems for easier selection.

The Drill Guide Block offers the possibility for parallel-straight and perpendicular attachment to the handle, adapting to anatomic requirements.



### Pre-Drilling Assembly

The Pre-Drilling Assembly consists of a Trocar, a Drill Sleeve and a Soft Tissue Protector which allows for pre-drilling and pin insertion without causing additional damage to the soft tissues. Different lengths enable you to choose the correct device for the soft tissue envelope. Dedicated Assemblies are available for 3mm, 4mm, 5mm and 6mm diameter pins.



**Note:**  
Do not tap on the trocar.

# Instruments

## Apex Instruments

Figure 1 illustrates manual Apex Pin Insertion using the Drill Brace, Drill Guide and Pre-Drilling Assemblies.



Fig .1

### Quick Release Apex Chuck

The Quick Release Apex Chuck is designed for fast and easy engagement of the Apex Pins and has a standard AO and a tri-flange connector. It is designed for insertion of Apex Pins by power tool.

### Combination Wrench/Pin Inserter

The Hoffmann II T-Wrench/Pin Inserter is used to insert 5mm & 6mm pins and tighten 7mm bolts. The Hoffmann II Compact Combination Wrench is used to insert 3mm & 4mm pins and tighten 5mm bolts.

### Note:

For the final seating in the second cortex the T-Wrench/Pin Inserter or the Drill Brace should be used.




# Technical Details

## Instructions for Instrument Use

### Drill Brace

The Drill Brace provides attachments for 3mm & 4mm pins on one end and 5mm & 6mm pins on the other end. For pin insertion, place the pin into the end correlating to the chosen pin diameter.

To access the different attachments for the pins, remove the handle and assemble it on the other end.





### Drill Guide Block

To assemble the Drill Guide Block, choose the correct block for your pin clamp. Set the Drill Guide Block in line or perpendicular and push it onto the handle aligning the laser-etched arrows. To release the block, push the button on the handle and pull it off.

### Quick Release Apex Chuck

To assemble the pin to the chuck, pull the sleeve backwards, as shown in the accompanying diagram and place the pin in the adapter. To secure the pin, push the sleeve back.

To release the pin from the adapter pull the sleeve backwards and remove the adapter from the pin.



# Technical Details

## Pin Insertion Guidelines

Among others, the Pin diameter influences axial frame rigidity. This is because the stiffness of the pin is a function of the forth power of the diameter<sup>5</sup>.

As a guideline one might use the following diameters:<sup>6</sup>



- Forearm: 3mm (distal) and 4mm (proximal) Apex Pins
- Humerus: 5mm Apex Pins, 4mm in distal fragments
- Femur and Pelvis: 5mm Apex Pins throughout the entire bone
- Tibia: 5mm Apex Pins
- Ankle: 5mm or 5mm transverse Apex Pins

The number of pins used in a frame construct depends on the patient condition and the indication. Increasing the number of pins will increase the frame rigidity. If available, correct position of the Apex Pins should always be verified by X-Ray examination.

When using a Self-Drilling/Self-Tapping Pin, turn the Drill Brace twice counter-clockwise to create a small notch for the pin. This helps prevent the pin from slipping on the cortex.

Afterwards, turn the Drill Brace clockwise for pin insertion.

Make a skin incision long enough and in the direction the skin will move during mobilization to avoid tension around the pin. This helps prevent irritation of the skin and may reduce the risk of infection.



5. Concepts in External Fixation, D. Seligson, 1982, Page 23 ff.

6. Bruce H. Ziran, Wade R. Smith, Jeff O. Anglen and Paul Tornetta, III; External Fixation: How to Make It Work; J Bone Joint Surg Am. 2007;89:1620-1632. Page 1627ff.

# Technical Details

## Pin Insertion Guidelines

When inserting the pin by power, using a low speed will limit the temperature increase, which can cause bone necrosis. Do not use excessive axial force. Figure 2 illustrates Apex Pin Insertion under power.

Insert the pins 90° to the long axis of the bone to reduce pull in and push out forces on the pins.





Fig. 2

However, there are additional factors such as fracture patterns and anatomical structures that need to be taken into consideration.

### Caution:

**Apex Pins are not intended for navigation procedure purposes.**




# Technical Details

## HA Apex Coated Pin Insertion Guidelines

Always use the Stryker Apex Predrilling Assemblies for pin insertion to avoid damage to the HA coating.

If the HA coating is damaged due to incorrect instrument usage, fixation properties may be compromised.



### Caution:

Do not wash or attempt to re-sterilize an unpacked HA coated Apex Pin. The coating may be damaged and the effect of the HA coating may be compromised.



Surgeons must always rely on their own clinical judgment when deciding which treatment and product to use with their patients.

# Ordering Information – Apex Pins

Product Number	Diameter mm Thread/Shaft	Total Length mm	Thread Length mm
<b>Self-Drilling/Self-Tapping</b>			
 5080-1-612	1.65/2.0	45	12
5080-1-620	1.65/2.0	45	20
5080-2-012	2.0	45	12
5080-2-020	2.0	45	20
5038-5-060	3.0	60	10
5038-1-080	3.0	80	10
5038-2-080	3.0	80	15
5038-5-080	3.0	80	20
5038-1-110	3.0	110	10
5038-2-110	3.0	110	25
5090-2-120*	3.0/4.0	120	20
5026-8-120	3.0/5.0	120	20
5023-1-090	4.0	90	10
5023-2-090	4.0	90	20
5023-3-090	4.0	90	30
5023-3-120	4.0	120	30
5023-5-120	4.0	120	35
5023-5-150	4.0	150	40
5023-6-150	4.0	150	50
5023-4-180	4.0	180	40
5023-6-180	4.0	180	50
5026-1-150	4.0/5.0	150	40
5018-3-120	5.0	120	30
5018-5-120	5.0	120	35
5018-5-150	5.0	150	40
5018-6-150	5.0	150	50
5018-3-180	5.0	180	35
5018-6-180	5.0	180	50
5018-8-180	5.0	180	60
5018-5-200	5.0	200	50
5018-6-200	5.0	200	60
5018-5-250	5.0	250	50
5018-7-250	5.0	250	70
5021-7-100	6.0	100	40
5021-7-150	6.0	150	50
5021-6-180	6.0	180	60
5021-8-200	6.0	200	70
5021-8-250	6.0	250	80
<b>Self-Drilling/Self-Tapping Titanium</b>			
 5039-2-110**	3.0	110	15
5016-5-111	5.0	120	35
5016-5-117	5.0	150	40
5016-5-118	5.0	150	50
5016-5-122	5.0	180	50

\* Available only NON-Sterile

\*\* Special order; available only NON-Sterile

**Note:**

Apex Pins are not intended for navigation procedure purposes.

# Ordering Information – Apex Pins

Product Number	Diameter mm Thread/Shaft	Total Length mm	Thread Length mm	
<b>Blunt/Self-Tapping</b>				
	5065-3-312	2.0	33	12
	5065-3-615	2.0	36	15
	5065-3-918	2.0	39	18
	5065-4-520	2.0	45	20
	5065-5-020	2.0	50	20
	5065-6-020	2.0	60	20
	5065-9-015	2.0	90	15
	5036-2-060	3.0	60	10
	5036-1-080	3.0	80	10
	5036-1-580	3.0	80	15
	5036-2-080	3.0	80	20
	5036-1-110	3.0	110	10
	5036-2-110	3.0	110	25
	5027-1-090	4.0	90	10
	5027-2-090	4.0	90	20
	5027-3-090	4.0	90	30
	5027-3-120	4.0	120	30
	5027-4-120	4.0	120	35
	5027-4-150	4.0	150	40
	5027-5-150	4.0	150	50
	5027-4-180	4.0	180	40
	5027-5-180	4.0	180	50
	5020-3-120	5.0	120	30
	5020-6-120	5.0	120	35
	5020-2-150	5.0	150	20
	5020-3-150	5.0	150	40
	5020-7-150	5.0	150	50
	5020-2-180	5.0	180	20
	5020-7-180	5.0	180	50
	5020-8-180	5.0	180	60
	5020-2-200	5.0	200	20
	5020-7-200	5.0	200	50
	5020-6-200	5.0	200	60
	5020-2-250	5.0	250	20
	5020-7-250	5.0	250	50
	5020-8-250	5.0	250	70
	5019-7-150	6.0	150	50
	5019-6-180	6.0	180	60
	5019-8-200	6.0	200	70
	5019-8-250	6.0	250	80

**Note:**

Apex Pins are not intended for navigation procedure purposes.

# Ordering Information – Apex Pins

	Product Number	Diameter mm Thread/Shaft	Total Length mm	Thread Length mm
<b>Cancellous</b>				
	5015-3-120	6.0/5.0	120	35
	5015-4-150	6.0/5.0	150	40
	5015-5-150	6.0/5.0	150	50
	5015-6-180	6.0/5.0	180	60
	5015-7-250	6.0/5.0	250	70
<b>Transfixing Pins</b>				
	5070-3-310*	2.0/1.5	33	10
	5070-3-810*	2.0/1.5	38	10
	5070-4-312*	2.0/1.5	43	12
	5070-4-812*	2.0/1.5	48	12
	5070-5-312*	2.0/1.5	53	12
	5070-5-815*	2.0/1.5	58	15
	5070-6-315*	2.0/1.5	63	15
	5070-7-820*	2.0/1.5	78	20
	5045-5-200**	3.0	200	N/A
	5030-3-200	5.0/4.0	200	35
	5030-4-200	5.0/4.0	200	40
	5030-5-200	5.0/4.0	200	50
	5030-3-250	5.0/4.0	250	35
	5030-4-250	5.0/4.0	250	40
	5030-5-250	5.0/4.0	250	50
	5030-6-250	5.0/4.0	250	60
	5030-4-300	5.0/4.0	300	40
	5030-5-300	5.0/4.0	300	50
	5030-7-300	5.0/4.0	300	70
	5050-5-250***	6.0/5.0	250	50
	5050-4-300*	6.0/5.0	300	40
	5050-5-300*	6.0/5.0	300	50

\* Available only NON-Sterile

\*\* Smooth Transfixing Pin Apex

\*\*\* Special order; available only NON-Sterile

**Note:**

Apex Pins are not intended for navigation procedure purposes.

# Ordering Information – HA Apex Pins

Stainless Steel REF	Diameter mm Thread/Shaft	Total Length mm	Thread Length mm
<b>HA Coated Self-Drilling/Self-Tapping</b>			
 5013-3-090S*	4.0	90	30
5013-2-120S*	4.0	120	20
5013-8-120S*	4.0	120	25
5013-3-120S*	4.0	120	30
5013-9-120S*	4.0	120	35
5013-2-150S*	4.0	150	20
5013-8-150S*	4.0	150	25
5013-3-150S*	4.0	150	30
5013-4-150S*	4.0	150	40
5017-9-120S	5.0	120	35
5017-2-150S	5.0	150	20
5017-8-150S	5.0	150	25
5017-3-150S	5.0	150	30
5017-9-150S	5.0	150	35
5017-4-150S	5.0	150	40
5017-5-150S	5.0	150	50
5017-6-150S	5.0	150	60
5017-3-180S	5.0	180	30
5017-4-180S	5.0	180	40
5017-5-180S	5.0	180	50
5017-6-200S*	5.0	200	60
5017-7-200S*	5.0	200	70
5014-2-120S*	6.0	120	30
5014-8-150S*	6.0	150	25
5014-3-150S*	6.0	150	30
5014-4-150S	6.0	150	40
5014-5-150S	6.0	150	50
5014-6-150S*	6.0	150	60
5014-3-180S	6.0	180	30
5014-4-180S	6.0	180	40
5014-5-180S	6.0	180	50
5014-6-180S*	6.0	180	60
5014-3-200S*	6.0	200	30
5014-4-200S	6.0	200	40
5014-5-200S	6.0	200	50
5014-6-200S	6.0	200	60
5014-7-250S*	6.0	250	70
5014-8-250S	6.0	250	80
5014-9-250S*	6.0	250	90

\* Special Order

**Note:**

Apex Pins are not intended for navigation procedure purposes.

All HA Apex Pins are single packed and sterile.

# Ordering Information – HA Apex Pins



Stainless Steel REF	Diameter mm Thread/Shaft	Total Length mm	Thread Length mm
<b>HA Coated Blunt</b>			
5008-2-150S	4.0	150	20
5008-8-150S	4.0	150	25
5008-3-150S	4.0	150	30
5008-4-150S	4.0	150	40
5009-2-200S	5.0	200	20
5009-8-200S	5.0	200	25
5009-3-200S	5.0	200	30
5009-9-200S	5.0	200	35
5009-4-200S	5.0	200	40
5009-5-200S	5.0	200	50
5009-6-200S*	5.0	200	60
5009-7-250S*	5.0	250	70
5009-8-250S*	5.0	250	80
5009-9-250S*	5.0	250	90
5010-8-200S*	6.0	200	25
5010-3-200S*	6.0	200	30
5010-4-200S	6.0	200	40
5010-5-200S	6.0	200	50
5010-6-200S*	6.0	200	60
5010-7-250S*	6.0	250	70
5010-8-250S*	6.0	250	80
5010-9-250S*	6.0	250	90

## Drill Bits for Apex Pins



5085-1-222	Drill Bit 2.2mm x 100mm for 3mm Pins
5085-2-032	Drill Bit 3.2mm x 200mm for 4mm Pins
5085-2-040	Drill Bit 4.0mm x 200mm for 5mm Pins
5085-2-045	Drill Bit 4.5mm x 200mm for 6mm Pins

\* Special Order


### Note:

Apex Pins are not intended for navigation procedure purposes.

# Ordering Information – Instruments

	Product Number	Description
	5057-0-300	Drill Brace Assembly
	5057-0-310	Handle for Drill Brace
	4920-9-030	7mm Wrench/5-6mm Pin Inserter
	4940-9-030	5mm Wrench/3-4mm Pin Inserter
	5057-1-003	Quick Release Apex Chuck with AO fitting, 3mm
	5057-1-004	Quick Release Apex Chuck with AO fitting, 4mm
	5057-1-005	Quick Release Apex Chuck with AO fitting, 5mm
	5057-1-006	Quick Release Apex Chuck with AO fitting, 6mm
	4922-9-050	Universal Chuck for 4, 5 & 6mm dia. Apex Pins, AO Coupling
	5057-1-110	Drill Guide Handle
	5057-1-115	Drill Guide Block, 5 hole, Hoffmann II MRI; Hoffmann 3, blue
	5057-1-116	Drill Guide Block, 10 hole, Hoffmann II MRI; Hoffmann 3, blue
	5057-1-117	Drill Guide Block, 4 hole, Hoffmann II Compact MRI, yellow
	5057-1-118	Drill Guide Block, Hoffmann II Compact Peri Articular Clamp, yellow
	5057-1-119	Drill Guide Block, 4 hole, Monotube Triax, blue, red
	5057-1-120	Drill Guide Block, 2 hole, Monotube Triax, yellow

# Ordering Information – Instruments

	Product Number	Description	Protection Length in MM	
	<b>APEX Instruments Pre-Drilling Assembly</b>			
	5057-3-100	Pre-Drilling Assembly, 3.0mm, short	33mm	
	5057-3-200	Pre-Drilling Assembly, 3.0mm, long	43mm	
	5057-4-000	Pre-Drilling Assembly, 4.0mm, extra short	35mm	
	5057-4-100	Pre-Drilling Assembly, 4.0mm, short	70mm	
	5057-4-200	Pre-Drilling Assembly, 4.0mm, long	100mm	
	5057-5-000	Pre-Drilling Assembly, 5.0mm, extra short	50mm	
	5057-5-100	Pre-Drilling Assembly, 5.0mm, short	73mm	
	5057-5-200	Pre-Drilling Assembly, 5.0mm, long	113mm	
	5057-6-000	Pre-Drilling Assembly, 6.0mm, extra short	60mm	
	5057-6-100	Pre-Drilling Assembly, 6.0mm, short	90mm	
	5057-6-200	Pre-Drilling Assembly, 6.0mm, long	120mm	
		4922-9-140	Tissue Protection Sleeve	
		4922-9-240	Trocar	
	5057-6-300	Pin Cutter, 4mm, 5mm & 6mm, Extension Handles		
<b>Disposable End Caps</b>				
	5047-1-030	Disposable End Caps 3mm, brown (15 pieces/pack)		
	5027-1-040	Disposable End Caps 4mm, white (15 pieces/pack)		
	5027-1-050	Disposable End Caps 5mm, blue (15 pieces/pack)		
	<b>APEX Storage Tray</b>			
	5057-9-913	Storage Tray Lid		
	5057-9-912	Storage Tray Upper Insert		
	5057-9-911	Storage Tray Lower Insert		
	5057-9-910	Storage Tray Base		

# Notes

# Notes

# Notes

## Reconstructive

---

Hips  
Knees  
Trauma & Extremities  
Joint Preservation  
Orthobiologics

## Medical & Surgical

---

Power Tools & Surgical Accessories  
Image Guided Navigation  
Endoscopy & Arthroscopy  
Integrated Communications  
Beds, Stretchers & EMS  
Sustainability Solutions

## Neurotechnology & Spine

---

Craniomaxillofacial  
Interventional Spine  
Neurosurgical, Spine & ENT  
Neurovascular  
Spinal Implants



### Manufactured by:

Stryker GmbH  
Bohnackerweg 1  
CH - 2545 Selzach  
Switzerland

[www.stryker.com](http://www.stryker.com)

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for Cleaning and Sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Apex, Hoffmann, Stryker, Tenxor, Triax. All other trademarks are trademarks of their respective owners or holders.

The products listed above are CE marked.

**Literature Number: 2012800**

**Content ID: APEX-ST-1, 08-2015**