


# Accolade® II

## Femoral Hip Stem

**Design rationale**



**Designed to fit more patients,  
designed to fit your approach.**



## Addressing modern demands with innovative technology.

The global THA population is evolving to include a younger,<sup>1</sup> more active<sup>2</sup> and more demanding<sup>3</sup> patient. Many femoral stem designs on the market today predate the emergence of this novel demand. Subsequently, an opportunity to enhance the conventional femoral stem design exists.

Conventional tapered wedge femoral stems have achieved popularity due to their simplicity and solid clinical results.<sup>4,5</sup> Despite these results, literature indicates that there are still unmet clinical needs.<sup>6,7,8</sup> Incidence of subsidence,<sup>6</sup> distal-only implant engagement<sup>7</sup> and periprosthetic fracture<sup>8</sup> suggest a clinical need for an improved implant fit for this novel patient population.


Stryker, along with key industry leaders, committed to developing an innovative femoral stem. This design would build from the sound principles of tapered wedge philosophy to meet the unique needs of the current patient population. At the heart of this development was a unique technology called Stryker Orthopaedic Modeling and Analytics, or SOMA. As a system that enables population-based design, SOMA has functionality with which to design, model and analyze novel orthopaedic devices.

Stryker utilized SOMA technology to design an innovative stem building upon the conventional tapered wedge femoral design, incorporating unique features to allow for an enhanced implant fit in today's patient population.<sup>9</sup> By establishing an increased canal fit and fill,<sup>9</sup> Accolade II has been shown to allow for improved stability,<sup>10</sup> decreased intraoperative femoral fractures,<sup>11</sup> as well as very good survivorship and functional outcomes,<sup>12,13</sup> ultimately leading to satisfied patients.<sup>13,14</sup>



## SOMA technology

Utilizing the proprietary SOMA technology, Stryker was able to complete one of the largest proximal femoral bone morphology studies ever undertaken.<sup>13</sup> An illustrated look at the process by which SOMA technology is employed in implant design is described below.




### CT acquisition

The SOMA database continues to acquire new CT scans and currently contains over 26,000 bones<sup>25</sup>




### Segmentation

Once acquired, all bones are segmented into inner and outer cortices




### Analysis

Using SOMA tools, bone morphology can be studied in a highly accurate and reproducible manner



### Design input

The results of these studies, such as the population canal flare index, can be utilized in implant design



### Validation

The resulting implant design can then be validated using SOMA fitting tools


## Three key SOMA-designed features

Bone morphology data allowed Stryker an unprecedented look at femoral anatomy and assisted in the design of an innovative femoral stem. The SOMA input\* was instrumental in the establishment of three key design features of Accolade II:



### Unique size-specific medial curvature

increasing proximal conformity to improve primary stability<sup>14,15</sup>



### Enhanced proximal-distal proportions<sup>16</sup>

shown to mimic canal anatomy to avoid distal-only engagement and achieve cortical fit<sup>8,16</sup>



### Optimized stem length\*\*

enables muscle-sparing approaches without sacrificing stability<sup>15,17</sup>

\*SOMA design of Accolade II based on 556 CT scans

\*\*Optimized stem length for a broad range of patients based on average femoral bone shape of 556 CT scans



## Unique size-specific medial curvature


Initial stability is critical to long-term implant performance.<sup>6</sup> Early subsidence and micromotion have been established as strong indicators for implant failure.<sup>6</sup> Initial stability may be increased by creating a higher conformity between the implant and the femoral cortices, leading to a larger area of contact.<sup>19</sup>

Analyzing the SOMA morphology study data,\* it was observed that a constant medial curvature may not allow for a conforming canal fit throughout varying femoral size.

This population-based input influenced Accolade II to incorporate the market's first unique size-specific medial curvature. This feature was designed to enable a more conforming proximal cortical fit,<sup>15</sup> which has been shown to allow for improved implant stability.<sup>19</sup>

Using femurs from the SOMA database,\* fit patterns of Accolade II can be compared to conventional tapered wedge designs. The three examples below illustrate how Accolade II achieves a more conforming canal fit throughout varying bone sizes.

### Small femur



Accolade II

Conventional tapered wedge

### Medium femur



Accolade II

Conventional tapered wedge

### Large femur



Accolade II


Conventional tapered wedge

\*SOMA design of Accolade II based on 556 CT scans


## Enhanced proximal-distal proportions

Distal-only engaged femoral stems can experience stress shielding,<sup>19</sup> and consequently may lead to elevated failure rates due to loosening and migration.<sup>9</sup> In order to better mimic the femoral anatomy and avoid distal-only engagement, a more anatomic implant growth rate is needed.

Utilizing the SOMA femoral morphology study,\* a more anthropometric proximal-distal stem growth rate was identified. This rate led to enhanced implant proportions,<sup>13</sup> as the distal geometry of Accolade II increases in size less than the proximal geometry. These proportions enable Accolade II to achieve a significantly better canal fit and fill<sup>9</sup> and have shown a decreased incidence of distal-only engagement.<sup>9</sup>



## Comparing implant fit



In the graph above, the proximal (P) and distal (D) measurements of a population of 556 femurs were plotted (black dots) against the corresponding stem diameters of Accolade II (gold) and a conventional tapered wedge design (gray). Accolade II achieved more fully conforming and proximal-only fit types compared to the conventional design, while subsequently **reducing distal-only fit by 14%.**<sup>16</sup>


\*SOMA design of Accolade II based on 556 CT scans

## Optimized stem length\*

Popularity of muscle-sparing approaches and bone-conserving fundamentals have led to a trend in shortening of stem length.<sup>20</sup> However, there is a complex relationship between stem length and implant stability.<sup>20</sup> Shortening stem length without geometry optimization has been shown to increase the potential for micromotion,<sup>20</sup> which is a strong indicator for implant failure.<sup>6</sup>

Accolade II utilized the SOMA database\*\* and stability analyses to establish an optimized length for each stem size, which not only accommodates muscle-sparing approaches,<sup>20</sup> but demonstrates improved initial stability.<sup>10</sup>

**“Simply shortening a standard tapered wedge design may reduce the primary stability.”<sup>17</sup>**



\*Optimized stem length for a broad range of patients based on average femoral bone shape of 556 CT scans

\*\*SOMA design of Accolade II based on 556 CT scans





## Designed to fit your approach


Muscle-sparing surgical approaches continue to gain popularity due to the potential patient benefits of faster recovery,<sup>21,22</sup> less pain<sup>21,22</sup> and greater satisfaction.<sup>23</sup>

Stryker's portfolio of muscle-sparing techniques features modern instrumentation and dynamic Medical Education programs to support the **Direct Anterior Approach** and the **Direct Superior Approach**.

The Direct Superior Approach was designed for surgeons who prefer the fundamentals and familiarity of the posterior approach, but seek to provide the next evolution in muscle-sparing THA surgical techniques for their patients.

Stryker has developed Stryker's Training Academy, a training platform to help surgeons manage the learning curve of a new surgical approach.

For access to Stryker's Training Academy, contact your local Stryker Sales Representative.



### DIRECT SUPERIOR APPROACH

**Toolkit** Return to menu

**Scenario: Femoral Neck Resection**  
 You are doing your first DSA case. The procedure has gone well so far. You have dissected and retracted the confluence of the piriformis and obturator internus tendons and performed the capsulotomy, but you run into difficulty during the osteotomy. While the hip may be dislocated, you still cannot adequately visualize if the resection has been made at the templated level.

**Patient profile:**  
 • 58-year-old female • BMI of 32

**What should you do? Make the best selection.**

**A Decision A**  
 Internally rotate the femur to 90 degrees.

**B Decision B**  
 Perform a subcapital osteotomy. Then replace the femoral retractors to visualize the remaining neck and make a second neck cut. Assess the combined neck resections to the pre-op templated neck resection height.

**C Decision C**  
 Reduce the head back into the acetabulum and make an in situ neck resection.

## Clinical performance

“Significantly better overall canal **fit**”<sup>9</sup>  
than conventional tapered wedge design


over  
**1,000,000**  
implanted worldwide<sup>24</sup>

**100%**  
**maintenance of bone  
mineral density**  
at the medial calcar at two years<sup>26</sup>

**100%**  
**aseptic survivorship**  
demonstrated in a midterm study<sup>12</sup>

**less than  
0.1 mm**  
**subsidence**  
observed in a two-year RSA study<sup>10</sup>

**7.5x**  
**fewer intraoperative  
fractures**  
observed compared to  
conventional tapered wedge<sup>11</sup>



## Accolade II Implant catalog numbers

Part number	Size	Neck angle
6720-0027	0	132°
6720-0127	1	
6720-0230	2	
6720-0330	3	
6720-0435	4	
6720-0535	5	
6720-0635	6	
6720-0737	7	
6720-0837	8	
6720-0937	9	
6720-1040	10	
6720-1140	11	

6721-0027	0	127°
6721-0127	1	
6721-0230	2	
6721-0330	3	
6721-0435	4	
6721-0535	5	
6721-0635	6	
6721-0737	7	
6721-0837	8	
6721-0937	9	
6721-1040	10	
6721-1140	11	

### References:

1. Health, United States, 2010: With Special Feature on Death and Dying. National Center for Health Statistics; 2011:17,60; DHHS Pub No. 2011-1232.
2. Kurtz SM, Lau E, Ong K, Zhao K, Kelly M, Bozic KJ. Future young patient demand for primary and revision joint replacement: national projections from 2010 to 2030. Clin Orthop Relat Res. 2009;467(10):2606-1612. doi:10.1007/s11999-009-0834-6
3. Haddad FS. A positive End to 2016. Bone Joint J. 2016;98-B(12):1569-1570. doi:10.1302/0301-620X.98B12.38083
4. Casper DS, Kim GK, Restrepo C, Parvizi J, Rothman RH. Primary total hip arthroplasty with an uncemented femoral component five- to nine-year results. J Arthroplasty. 2011;26(6):838-841. doi:10.1016/j.arth.2011.02.010
5. McLaughlin JR, Lee KR. Total hip arthroplasty with an uncemented tapered femoral component. J Bone Joint Surg Am. 2008;90(6):1290-1296. doi:10.2106/JBJS.G.00771
6. van der Voort P, Pijls BG, Nieuwenhuijse MJ, et al. Early subsidence of shape-closed hip arthroplasty stems is associated with late revision. A systematic review and meta-analysis of 24 RSA studies and 56 survival studies. Acta Orthop. 2015;86(5):575-585. doi:10.3109/17453674.2015.1043832
7. Cooper HJ, Jacob AP, Rodriguez JA. Distal fixation of proximally coated tapered stems may predispose to a failure of osteointegration. J Arthroplasty. 2011;26(6 Suppl):78-83. doi:10.1016/j.arth.2011.04.003
8. Abdel MP, Watts CD, Houdek MT, Lewallen DG, Berry DJ. Epidemiology of periprosthetic fracture of the femur in 32 644 primary total hip arthroplasties: a 40-year experience. Bone Joint J. 2016;98-B(4):461-467. doi:10.1302/0301-620X.98B4.37201
9. Issa K, Pivec R, Wuestemann T, Tatevossian T, Nevelos J, Mont MA. Radiographic fit and fill analysis of a new second-generation proximally coated cementless stem compared to its predicate design. J Arthroplasty. 2013; 29(1):192-198. doi:10.1016/j.arth.2013.04.029
10. Collopy D. A prospective single-centre case series using Roentgen stereophotogrammetric analysis (RSA) to evaluate stem micromotion of the Stryker Accolade II cementless hip stem to two years post-surgery. Bone Joint J:Orthop Proc. 99-B(Suppl\_3):77.
11. Fleischman AN, Schubert MM, Restrepo C, Chen AF, Rothman RH. Reduced incidence of intraoperative femur fracture with a second-generation tapered wedge stem. J Arthroplasty. 2017;32(11):3457-3461. doi:10.1016/j.arth.2017.06.018
12. Kolisek FR, Jaggard CE, Milto AJ, et al. Mid-term follow up of newer-generation morphometric wedge stems for total hip Arthroplasty (THA). Surg Technol Int. 2020;36:399-403. Based on 186 THA patients with average follow-up of 65.7 months.
13. Kolisek FR, Chughtai M, Mistry JB, et al. Outcomes of second-generation tapered wedge design. Surg Technol Int. 2016; 28:275-279.
14. Stryker internal presentation. Global bone morphology study. 2016.
15. Narzikul A, Restrepo C, Chen AF, Parvizi J. Alteration in geometry of femoral stem results in better fit and fill: comparison of Accolade® I vs Accolade® II. Presented at: Eastern Orthopaedic Association 47th Annual Meeting; October 19-22, 2016; New Orleans, LA.
16. Wuestemann T, Bastian A, Parvizi J, Nessler J, Kolisek F, Nevelos J. A novel tapered hip stem design optimized for femoral fit in a wide array of bone types. Presented at: 12th EFORT Annual Congress; June 1-4, 2011; Copenhagen, Denmark.
17. Boucher F, Wuestemann T, Parvizi J, Nevelos J. Preclinical finite element analysis of a novel tapered wedge stem optimized for stability at a reduced length. Presented at: 12th EFORT Annual Congress; June 1-4, 2011; Copenhagen, Denmark.
18. Lovell T, Hozack W, Kreuzer S, et al. Influence of stem length on the insertion path in THR. Presented at: 55th Annual Meeting of the Orthopaedic Research Society (ORS); February 22-25, 2009; Las Vegas, NV.
19. Khanuja HS, Vakil JJ, Goddard MS, Mont MA. Cementless femoral fixation in total hip arthroplasty. J Bone Joint Surg Am. 2011;93(5):500-509. doi:10.2106/JBJS.J.00774
20. Faizan A, Wuestemann T, Nevelos J, Bastian AC, Collopy D. Development and verification of a cementless novel tapered wedge stem for total hip arthroplasty. J Arthroplasty. 2015;30(2):235-240. doi:10.1016/j.arth.2014.09.023
21. Restrepo C, Parvizi J, Eslam Pour AE, Hozack WJ. Prospective randomized study of two surgical approaches for total hip arthroplasty. J Arthroplasty. 2010;25(5):671-679.e1. doi:10.1016/j.arth.2010.02.002
22. Vail T, Mariani EM, Bourne MH, Berger RA, Meneghini RM. Approaches in primary total hip arthroplasty. J Bone Joint Surg Am. 2009;91(Suppl 5):10-12. doi:10.2106/JBJS.I.00343
23. Zawadsky MW, Paulus MC, Murray PJ, Johansen MA. Early outcome comparison between the direct anterior approach and the mini-incision posterior approach for primary total hip arthroplasty: 150 consecutive cases. J Arthroplasty. 2014; 29(6):1256-1260. doi:10.1016/j.arth.2013.11.013
24. Data on file. Stryker sales data, 2020.
25. Stryker internal memo
26. Nam D, Salih R, Nahhas CR, Barrack RL, Nunley RM. Is a modular dual mobility acetabulum a viable option for the young, active total hip arthroplasty patient? Bone Joint J. 2019;101-B(4):365-371. doi:10.1302/0301-620X.101B4.BJJ-2018-0834.R1 Based on 43 patients with a minimum follow-up of 2 years

## Joint Replacement

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker's product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any of Stryker's products. The products depicted are CE marked according to the Medical Device Regulation 2017/745 or the Medical Device Directive 93/42/EEC. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your sales representative if you have questions about the availability of products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Accolade, Stryker, SOMA. All other trademarks are trademarks of their respective owners or holders.

ACCII-BRO-2\_Rev-1\_24634

Copyright © 2020 Stryker



### Manufactured by:

Howmedica Osteonics Corp.  
325 Corporate Drive  
Mahwah, NJ 07430, USA  
A subsidiary of Stryker Corp.  
**stryker.com**