
*s*tryker

## **ANATO**®

**Femoral Hip System** 

-Surgical Protocol



#### **INDICATIONS**

The indications for use of the total hip replacement prostheses include:

- Noninflammatory degenerative joint disease, including osteoarthritis and avascular necrosis;
- Rheumatoid arthritis;
- Correction of functional deformity;
- Revision procedures where other treatments or devices have failed; and,
- Nonunions, femoral neck fractures, and trochanteric fractures of the proximal femur with head involvement that are unmanageable using other techniques.

Additional indication specific to use of Anato Stems are with compatible Howmedica Osteonics Constrained Liners:

 When the stem is to be used with compatible Howmedica Osteonics Constrained Liners, the device is intended for use in primary or revision patients at high risk of hip dislocation due to a history of prior dislocation, bone loss, soft tissue laxity, neuromuscular disease, or intra-operative instability.

The Anato Hip Stem is intended for cementless use only and is intended for total and hemi-arthroplasty procedures.

#### **CONTRAINDICATIONS**

- Active infection or suspected latent infection in or about the hip joint;
- Bone stock that is inadequate for support or fixation of the prosthesis;
- Skeletal immaturity; and,
- Any mental or neuromuscular disorder that would create an unacceptable risk of prosthesis instability, prosthesis fixation failure, or complications in post-operative care.

#### **WARNINGS AND PRECAUTIONS**

See package insert for warnings, precautions, adverse effects and other essential product information.

Before using Anato Stem instrumentation, verify:

- Instruments have been properly disassembled prior to cleaning and sterilization;
- Instruments have been properly assembled post-sterilization;
- Instruments have maintained design integrity; and,
- Proper size configurations are available.

For Instructions for Cleaning, Sterilization, Inspection, and Maintenance of Orthopaedic Medical Devices, refer to LSTPI-B and L24002000.

#### **TABLE OF CONTENTS**

Pre-operative Planning
Femoral Neck Resection
Preparing the Femoral Canal
Broaching6-7
Trial Reduction
Implanting the Stem
Final Reduction
Implant Information
Appendix
Instrumentation



#### **Femoral Hip System**

The Anato Hip Stem is intended for cementless, press-fit application. The proximal region of the stem is coated with PureFix HA over a commercially pure titanium plasma spray substrate. The Anato femoral hip system is suitable for various surgical approaches, including direct anterior.

The total system includes:

- Left and right monolithic femoral hip implants
- Neutral and anteverted neck options. There is a 7 degree difference between the two options
- 8 body sizes ranging from size 1 to size 8
- One anatomic stem offset of 130 degrees

The stem is designed for use with Stryker V40 femoral heads and their compatible acetabular components.

This publication sets forth detailed recommended procedures for using Stryker Orthopaedics devices and instruments. It offers guidance that a surgeon should need, but, as with any such technical guide, each surgeon must consider the particular needs of each patient and make appropriate adjustments when and as required.

#### PRE-OPERATIVE PLANNING








Figure 1

Pre-operative planning aids in

the determination of probable implant style and size. The preoperative planning process should take qualitative and quantitative factors (including patient bone quality, density, and morphology) into consideration in order to evaluate and select the appropriate instrument/implant system for the patient.

When templating the acetabulum, place an acetabulum template over the area on the X-ray. Be sure that the cup is well centered within the acetabulum and the size fills between the tear drop and the superior rim. After templating the acetabulum for size, mark the center of rotation represented by the purple dot (Figure 1).

Optimal femoral stem fit, prosthetic neck length, angle and version can be evaluated with the use of pre-operative X-ray analysis.

The following parameters should be determined using an A/P radiograph: stem size, femoral offset, leg length, neck angle, and center of rotation.

The Anato stem provides offset across all sizes through a 130° neck angle. Choose the Anato template for which the stem size achieves adequate metaphyseal fill and recreates the desired leg length and offset (Figure 2). This engagement can be seen when the end of the lateral flare is at the transition between the metaphyses and diaphyses (see arrow in Figure 2). The proximal end of the coated region should then be at the level of the piriformis fossa. The template has markings that indicate the center of the femoral head for

Figure 3

a range of head offset options. The predicted change in leg length and offset is determined by the relative positioning of the center of rotation markings on the femoral and acetabular components. For example, if a given femoral component center of rotation marking is superior to the center of rotation marking of the acetabular component, leg lengthening is predicted. The desired change in leg length is determined by the radiographic leg length inequality that was previously determined with an A/P bilateral film. If 8mm of leg lengthening is required in order to equalize the leg lengths, the center of rotation marking of the femoral component must be 8mm above the center of rotation of the acetabular component. Once the final estimated stem size and position is determined, the neck resection level should be noted (Figure 3). This will be used as a reference during intra-operative neck resection.



**ANATO** Acetate **Templates** Lit# LTEM107




Acetate templates include 20% magnification.

#### **FEMORAL NECK RESECTION**







A proper neck resection level directly affects the final placement and fit of the femoral stem. As the visual landmarks for surgical approaches may vary, it is important to note that different landmarks may be used to aid in locating the appropriate resection. (Figure 4). The height of the cut can be adjusted in line with the templating activity in order to reconstruct the correct patient anatomy.

Illustrated is a view of the proximal femur with the anterior landmarks highlighted (Figure 4). Outlined in blue is the intertrochanteric line which represents the origin of the quadriceps muscle and the insertion of the anterior capsule. The second landmark, the saddle (green line), is a shape created by the outline of the superior aspect of the head-neck-greater trochanter. This helps define the starting point for the femoral neck osteotomy.

#### **Single Osteotomy**

The neck osteotomy should start at the saddle proximal to the intertrochanteric line and extend distally at approximately a 45 degree angle with respect to the femoral shaft (Figure 5). If the intertrochanteric line cannot be approximated, an estimation of the femoral axis will also aid in approximating the direction of the cut.

#### **Dual Osteotomy (Optional)**

Illustrated is a view of the proximal femur highlighting the landmarks that are typically visible from the posterior approach (Figure 6). The point represents the insertion area of the piriformis tendon. This neck resection consists of two separate cuts. The first is a vertical cut to the depth of the insertion of the piriformis tendon. The second cut is approximately 60 degrees from the femoral axis to the medial resection cut.



#### Note

The medial resection point will serve as a reference for depth and seating of the broaches and stems in the subsequent steps and operations.

#### PREPARING THE FEMORAL CANAL

To help ensure proper final orientation of the stem, lateral preparation may be required. Retraction of the gluteus medius and removal of a portion of the superior lateral neck and/or piriformis fossa will permit true axial introduction of the instruments and implant. A box osteotome is useful in removing bone from this area. Connect the box osteotome to the broach handle. Use a mallet to strike the end of the broach handle with careful controlled blows to remove bone of the superior femoral neck and open the medullary canal.





MIS Modular 8mm Hollow Box Chisel 4849-8-108



Offset Broach Handle 1020-1460



#### PREPARING THE FEMORAL CANAL CONTINUED



#### **Distal Reaming**

The Anato hip stem design encourages metaphyseal load transfer. This design enables engagement in the metaphysis and the prevention of distal fixation. If the pre-operative planning anticipates a potential contact between the stem and the diaphyseal cortex then it is recommended to ream with flexible reamers up to the diameter corresponding to the templated implant size in Table 1 below. This decision may also be made at the time of broaching.

#### **TABLE 1: RECOMMENDED REAMING DIAMETER**

STEM SIZE	DISTAL STEM DIAMETER (MM)	RECOMMENDED REAMING DIAMETER (MM)
1	7.6	10.0
2	7.9	10.0
3	8.4	10.0
4	9.4	11.0
5	10.4	12.0
6	11.7	13.0
7	13.5	14.0
8	15.2	16.0

Figure 8



#### Note

The BixCut IM fixed-head system with modified trinkle fitting, can be used to prepare the femur to receive an Anato hip stem. Be sure to use the 0227 BixCut family of reamers. (See instrument list in the Appendix for full list of part numbers).



Instrument

BixCut Reamer 6mm - 16mm 0227-XXXX



#### Note

It is recommended that the starting reamer diameter be 2mm smaller than the final recommended reaming diameter as described in column three of Table 1 above.

#### **Guide Wire Information:**

For BixCut Reamers sizes 6 & 7, the compatible guide wire is 2.5mm in diameter (part #1806-0084S).

All other BixCut Reamers are compatible with guide wires with a 3.0mm diameter (part #1806-0085S).

#### **BROACHING**



Assemble the broach to the broach handle. An audible click will be heard when the broach is secured onto the broach handle.

Begin with the size 1 broach and progressively increase the broach size until the appropriate amount of size is achieved (Figure 9).



Note

Additional bone preparation may be performed with standard surgical instruments, such as currettes, rongeurs, etc.



Note

Apply a lateral force to the broach during impaction to avoid a varus position. The broach will determine the size and the position of the final implant.



Note

The lateral level of the broach may sit below the level of the neck resection laterally.

Figure 9

Appropriate fill will be denoted by a change in pitch, tactile resistance, or when the broach ceases to advance. Verify a secure fit by ensuring that the broach has axial and rotational stability. With proper cortical contact, the broach should not twist or turn relative to the femur. If there is movement, a larger size broach may be needed, or distal reaming may be necessary as the broach may be contacting distally in diaphyseal bone.

Generally, if the broach sinks below the level of the neck resection, advance to the next larger broach. Relying on the neck resection height alone for final seating height may lead to improper sizing and inadequate stem fixation. The final broach should seat firmly in the metaphysis.

Upon reaching the final size and depth of the broach, detach the broach handle from the broach, leaving the broach fully seated in the femoral canal.



Note

If the broach is unstable either before or after trial reductions the use of a larger broach is recommended. If it is suspected that the broach is hanging up distally, distal reaming may be required.



Tip

The appropriate size is achieved when the broach stops advancing and the shoulder seats at approximately the same level as the trochanteric fossa. The last tooth of the broach corresponds to the upper limit of the coating on the stem as seen in Figure 10.




**Instruments** 

Offset Broach Handle 1020-1460



Anato Broach Size 1-8 Left/Right 4845-2-9XX 4845-7-9XX

#### **BROACHING** CONTINUED



#### **OPTIONAL STEP**

Take an intra-operative x-ray if needed to evaluate alignment and fit of the selected implant.

#### **Optional Instruments:**


The existing ABG II Monolithic hip stem broach sizes 6-8 are compatible with this system, since they provide the same proximal fit.

PART NUMBER	DESCRIPTION
4845-2-956	ABGII MIS Broach Size 6 Right
4845-2-957	ABGII MIS Broach Size 7 Right
4845-2-958	ABGII MIS Broach Size 8 Right
4845-2-966	ABGII MIS Broach Size 6 Left
4845-2-967	ABGII MIS Broach Size 7 Left
4845-2-968	ABGII MIS Broach Size 8 Left

Note that the ABG II Broaches Size 6-8 may remove additional bone distal to the implant in patients with tighter distal femoral canals.

Please refer to the appendix for additional broach handle information.

#### TRIAL REDUCTION



### Selection of neck trial (anteverted or neutral) and assembly of neck trial onto broach

The trial assembly, which consists of the broach, trial neck, and trial head (Figure 11), is used to provide an evaluation of the hip mechanics during a trial reduction. Before the selection and implantation of the final component, modifications to the pre-operative plan, in terms of neck length and/or head diameter, can be made. With the broach left in the femur, the surgeon can adjust the biomechanics through version options and head offsets. These can be balanced to achieve good internal/external rotational mobility in extension and good stability of the hip in flexion/internal rotation/adduction and in extension/external rotation.

Based on the chart below, select a neck trial which has the same version as the planned monolithic implant option.

#### **TABLE 2: NECK TRIAL SELECTION GUIDE**

COLOR	VERSION RIGHT OR LEFT	
Silver	Anteverted (7 degrees)	Right/Left
Black	Neutral	Right
Gold	Neutral	Left



#### Note

The same neck trials can be used for all implant sizes because the implant neck lengths are built into the broaches.



Instruments

Anato Left Neutral Neck Trial 4845-7-970



Anato Right Neutral Neck Trial 4845-7-971



Anato Left/Right Anteverted Neck Trial 4845-7-973



#### TRIAL REDUCTION CONTINUED



Assemble the neck trial onto the broach. Next, assemble a V40 head trial onto the neck trial (Figure 12).

Femoral heads come in multiple offsets and are different for each femoral head implant material and diameter (see Table 3). For this reason, final head material should be chosen prior to trial reduction. Offsets add or subtract from the base neck length of the implant and help to achieve the desired leg length and offset.

After the trial has been assembled to the broach, perform a trial reduction of the hip to evaluate the biomechanics. Upon confirmation of the selected components, remove the trial head and trial neck, and reassemble the broach handle to the broach and remove it from the femoral canal. The final broach size determines the correct implant size.



Ensure that the broach post and the V40 trunnion are clear of debris prior to assembling the neck and head trials.

**TABLE 3: HEAD COMPATIBILITY GUIDE** 

HEAD	HEAD SIZE	HEAD OFFSETS
	22	+0, +3, +8
	26	-3, +0, +4, +8
	28	-4, +0, +4, +8
CoCr V40	32	-4, +0, +4, +8
	36	-5, +0, +5
	40	-4, +0, +4, +8
	44	-4, +0, +4, +8
	28	-2.7, +0, +4
Alumina V40	32	-4, +0, +4
	36	-5, +0, +5
Alumina C-Taper	28	-2.5, +0, +5
(when used with C-Taper Adaptor Sleeve-catalog	32	-2.5, +0, +5
#17-0000E)	36	-5, +0, +5
	28	-4, -2.7, +0, +4
delta BIOLOX V40	32	-4, +0, +4
	36	-5, -2.5, +0, +2.5, +5, +7.5
delta BIOLOX C-Taper	28	-2.5, +0, +2.5, +5
(when used with C-Taper Sleeve	32	-2.5, +0, +2.5, +5
- catalog #17-0000E)	36	-5, -2.5, +0, +2.5, +5, +7.5
	28	-2.5, +0, +4
delta BIOLOX Universal Taper	32	-2.5, +0, +4
(when used with Universal Taper Sleeve -	36	-2.5, +0, +4
catalog #6519-T-XXXX)	40	-2.5, +0, +4
	44	-2.5, +0, +4

- \*+0mm is the maximum offset allowed for the size 1 stem
- \*+6mm is the maximum offset allowed for the size 2 stem
- \*+8mm is the maximum offset allowed for the size 3-8 stems



#### Note

For additional acetabular and head compatibility information please see Table 8 in the Appendix.

For additional instruction on acetabular and head implantation, please see the respective surgical technique for the compatible products referenced above.



#### Note

Head trials with an "R" suffix are made from a radiopaque material, making them visible on X-ray. Instrument 6264-x-xxxR V40 head trial



Instruments

V40 Head Trial 6264-x-xxxR



#### **IMPLANTING THE STEM**

The stem can be inserted by hand until approximately 1 cm of the coated region lies above the final broach position. Final implant seating is achieved by gently impacting the stem using the stem impactor (Figure 13). The use of the bullet tipped femoral impactor is intended to allow for rotation of the stem to follow the shape of the femoral canal.





#### Note

The surgeon should NOT attempt to continue impacting the femoral component if visual and auditory clues indicate that the stem is firmly seated in the canal. These clues, rather than the broach seating level, should be used to determine the final seating height of the implant. Continued aggressive impaction could lead to femoral fracture. In the event that dense bone is encountered intra-operatively and compounding anatomical factors are present, the seating of the implant may not be consistent with the level of the broach due to the viscoelastic nature of the femoral bone.



Instruments

Stem Impactor 4842-1-000



#### **IMPLANTING THE STEM CONTINUED**



Figure 14




Prior to implanting the final stem and femoral head, check the implants for potential damage. Ensure that the final stem matches the last broach used. Use the markings on the face of the trunnion to verify the stem size and taper type.



#### **TABLE 4: IMPLANT PART NUMBER REFERENCE GUIDE**

CATALOG NUMBER	GUIDE
	X indicates LEFT (2) or RIGHT (1) Implant
4845-7-XYZ	Y indicates Version –Anteverted (0) or Neutral (1)
	Z indicates Size of Implant (1-8)

#### FINAL REDUCTION



Prior to final head assembly, neck length/head offset may be re-evaluated using a V40 head trial. Place the head trial on to the stem neck taper and reduce the hip. Evaluate leg length and soft tissue tension. Remove the head trial and dry the implant trunnion with a laparotomy sponge or sterile towel.

Select the appropriate corresponding V40 femoral head or sleeve and place it onto the dry trunnion of the femoral stem with a slight twist. Impact the head with moderate impactions using the head impactor (Figure 15).

Verify the head is secure on the trunnion after head impaction by applying traction to the head confirming stability on the trunnion.

The use of a V40 universal adaptor sleeve is necessary when using a BIOLOX delta universal taper ceramic femoral head with an Anato stem. After completing the trialing process, intra-operatively assemble the adaptor sleeve to the femoral stem manually. The universal adaptor sleeve must be fully seated on the stem trunnion before the head is assembled.



#### **Note**

Ensure that the V40 taper is clean prior to assembling the head trial or head implant



#### Warnings

In no instance should any attempt be made to preassemble the adaptor sleeve inside the BIOLOX delta universal ceramic head. Intra-operatively assemble the BIOLOX delta universal taper ceramic head onto the sleeved femoral stem and set with moderate blows using the Stem Head Impactor (1104-1000). Care must be taken to avoid excessive impact forces when assembling the ceramic head to the sleeved femoral component.

The wound must be fully lavaged prior to wound closure in order to ensure the removal of any debris that was potentially generated during surgery.


Reduce the femoral head into the acetabular cup and re-check for range of motion and soft tissue tension. The surgical site is then closed according to surgeon preference.

### TABLE 5: UNIVERSAL ADAPTOR SLEEVES – TITANIUM

CATALOG NUMBER	OFFSET (MM)	TAPER
6519-T-025	-2.5	V40
6519-T-100	+0	V40
6519-T-204	+4	V40



### FINAL REDUCTION CONTINUED



#### **OPTIONAL STEP: Extraction**

When using the stem extractor, secure it below the neck so that it is oriented beneath the taper as shown (Figure 16). Mallet blows can be used on the underside of the strike plate to remove the stem. If additional force is required, the use of the McReynolds slap hammer may be used to aid in removing the stem. This instrument can be used in conjunction with the McReynolds slap hammer by connecting it to 6869-1-000 (in conjunction with 6869-2-000 and 6869-3-000).

#### **Optional Instruments:**

PART NUMBER	DESCRIPTION
4845-7-530	Stem extractor
6869-1-000	McReynolds Driver
6869-2-000	McReynolds Driver
6869-3-000	McReynolds Driver



#### Warnings

The use of this device may cause damage to the implant and reduce its service life. As such, if the stem extractor is used to remove an implant, the implant must not be reused.






#### **IMPLANT INFORMATION**



Stem Length is measured from the top of the lateral shoulder to the distal tip.

#### **TABLE 6: ANATO HIP STEM OFFSET TABLE**

				А	В	С	D
CATALOG NUMBER		STEM DESCRIPTION	SIZE	NECK ANGLE DEGREE SIZE	NECK LENGTH (MM)	STEM LENGTH (MM)	FEMORAL OFFSET WITH +0MM HEAD
	4845-7-101	Size 1 Right Anteverted	1		28	99	37
	4845-7-102	Size 2 Right Anteverted	2		28	99	39
	4845-7-103	Size 3 Right Anteverted	3		32	110	42
	4845-7-104	Size 4 Right Anteverted	4		32	115	44
	4845-7-105	Size 5 Right Anteverted	5		34	120	44
	4845-7-106	Size 6 Right Anteverted	6		34	120	47
ANTEVERTED	4845-7-107	Size 7 Right Anteverted	7		36	120	50
ER	4845-7-108	Size 8 Right Anteverted	8	130°	38	120	53
Ä	4845-7-201	Size 1 Left Anteverted	1		28	99	37
A	4845-7-202	Size 2 Left Anteverted	2		28	99	39
	4845-7-203	Size 3 Left Anteverted	3		32	110	42
	4845-7-204	Size 4 Left Anteverted	4		32	115	44
	4845-7-205	Size 5 Left Anteverted	5		34	120	44
	4845-7-206	Size 6 Left Anteverted	6		34	120	47
	4845-7-207	Size 7 Left Anteverted	7		36	120	50
	4845-7-208	Size 8 Left Anteverted	8		38	120	53
	4845-7-111	Size 1 Right Neutral	1		28	99	37
	4845-7-112	Size 2 Right Neutral	2		28	99	39
	4845-7-113	Size 3 Right Neutral	3		32	110	42
	4845-7-114	Size 4 Right Neutral	4		32	115	44
	4845-7-115	Size 5 Right Neutral	5		34	120	44
	4845-7-116	Size 6 Right Neutral	6		34	120	47
_	4845-7-117	Size 7 Right Neutral	7		36	120	50
RA	4845-7-118	Size 8 Right Neutral	8	130°	38	120	53
NEUTRAL	4845-7-211	Size 1 Left Neutral	1		28	99	37
Z	4845-7-212	Size 2 Left Neutral	2		28	99	39
	4845-7-213	Size 3 Left Neutral	3		32	110	42
	4845-7-214	Size 4 Left Neutral	4		32	115	44
	4845-7-215	Size 5 Left Neutral	5		34	120	44
	4845-7-216	Size 6 Left Neutral	6		34	120	47
	4845-7-217	Size 7 Left Neutral	7		36	120	50
	4845-7-218	Size 8 Left Neutral	8		38	120	53





#### **APPENDIX**

The Anato hip stem utilizes the same broach handle connection as the Accolade TMZF (LASST), Accolade II (LSP76) and ABG II Monolithic hip stem (LSP77). This is also the same broach handle connection as seen in the direct anterior instrument set (DAA-SS-1).

#### **TABLE 7: BROACH HANDLE OPTIONS**

PART DESCRIPTION	PART NUMBER
Offset Broach Handle	1020-1460
Straight Broach Handle	1440-1460

### TABLE 8: ANATO HIP STEM COMPATIBLE FEMORAL HEADS

HEAD	HEAD DIAMETER (MM)	OFFSET (MM)
V40 BIOLOX delta	28	-4.0, -2.7, +0, +4.0
V40 BIOLOX delta	32	-4.0, +0, +4.0
V40 BIOLOX delta	36	-5, -2.5, +0, +2.5, +5, +7.5
Universal Taper, BIOLOX delta Head	28, 32, 36, 40, 44	
Universal Taper Sleeves - V40 stems		-2.5, +0, +4
V40 CoCr	22	+0, +3, +8
V40 CoCr	26	-3, +0, +4, +8
V40 CoCr	28	-4, +0, +4, +8
V40 CoCr	32	-4, +0, +4, +8
V40 LFIT CoCr	22	+0, +3, +8
V40 LFIT CoCr	26	-3, +0, +4, +8
V40 LFIT CoCr	28	-4, +0, +4, +8
V40 LFIT CoCr	32	-4, +0, +4, +8
V40 LFIT CoCr	36	-5, +0, +5
V40 LFIT CoCr	40	-4, +0, +4, +8
V40 LFIT CoCr	44	-4, +0, +4, +8
Alumina V40 Head	28	-2.7, +0, +4
Alumina V40 Head	32	-4, +0, +4
Alumina V40 Head	36	-5, +0, +5
C-Taper Alumina	28, 32	-2.5, +0, +5
C-Taper Alumina	36	-5, +0, +5
C-Taper Delta	28, 32	-2.5, +0, +2.5, +5
C-Taper Delta	36	-5, -2.5, +0, +2.5, +5, +7.5
Unitrax V40 Head	38, 40-56, 58, 61	+0
Unitrax V40 Monolithic Adapter		-4, +0, +4, +8
V40/C-taper Adapter Sleeve		

- +0mm is the maximum offset allowed for the size 1 stem
- +6mm is the maximum offset allowed for the size 2 stem
- +8mm is the maximum offset allowed for the size 3-8 stems



#### **INSTRUMENTATION**

#### 4845-7-602 Anato Top Tray General Instruments

— 4845-7-530 Stem Extractor

> 4849-8-108 MIS Modular 8mm Hollow Box Chisel

4845-7-970


Neutral Neck Trial – Left

4845-7-971

Neutral Neck Trial - Right

4845-7-973

Neck Trial - Anteverted



4845-5-520 Head / Neck Impactor

#### Optional

4845-5-525 Curved Head/Neck Impactor

4842-1-000 Stem Impactor 6264-8-028R

28mm -4mm V40 Head Trial

6264-8-928R

28mm -2.7mm V40 Head Trial

6264-8-128R

28mm +0mm V40 Head Trial

6264-8-228R

28mm +4mm V40 Head Trial

6264-8-328R

28mm +8mm V40 Head Trial

6264-8-428R

28mm +12mm V40 Head Trial

6264-8-032R

32mm -4mm V40 Head Trial

6264-8-632R

32mm -2.5mm V40 Head Trial

6264-8-132R

32mm +0mm V40 Head Trial

6264-8-232E

32mm +4mm V40 Head Trial

6264-8-332R

32mm +8mm V40 Head Trial

6264-8-432R

32mm +12mm V40 Head Trial

6264-8-036R

36mm -5mm V40 Head Trial

6264-8-436R

36mm -2.5mm V40 Head Trial

6264-8-136R

36mm +0mm V40 Head Trial

6264-8-536R

36mm +2.5mm V40 Head Trial

6264-3-236R

36mm +4mm V40 Head Trial

6264-8-236R

36mm +5mm V40 Head Trial

6264-8-736R

36mm +7.5mm V40 Head Trial



#### **INSTRUMENTATION**

#### 4845-7-603 Anato Bottom Tray BixCut Reamers

0227-6080 BixCut Reamer 8mm

DIXOUT TOUTHOR OTHER

0227-6090 BixCut Reamer 9mm

0227-6100

BixCut Reamer 10mm

0227-6110

BixCut Reamer 11mm

0227-8120

BixCut Reamer 12mm

0227-8130

BixCut Reamer 13mm

0227-8140

BixCut Reamer 14mm

0227-8150

BixCut Reamer 15mm

0227-8160


BixCut Reamer 16mm

0227-5060

BixCut Reamer 6mm

0227-5070

BixCut Reamer 7mm



1806-0080S 3.0mm x 800mm Ball Tip Guide Wire 1806-0083S 2.5mm x 800mm Ball Tip Guide Wire



#### INSTRUMENTATION

4845-7-601 **ANATO Broach Tray** 4845-7-600 **Stryker Case** 

4845-2-961

ANATO Broach Size 1 Left

4845-2-962

ANATO Broach Size 2 Left

4845-2-963

ANATO Broach Size 3 Left

4845-2-964

ANATO Broach Size 4 Left

4845-2-965

ANATO Broach Size 5 Left

4845-7-966

ANATO Broach Size 6 Left

4845-7-967

ANATO Broach Size 7 Left

4845-7-968

ANATO Broach Size 8 Left

The longer broaches mentioned in step 4, broaching, can also fit into this tray in the place of the comparable size and side broach.

4845-2-956 ABGII MIS Broach Size 6 Right 4845-2-957 ABGII MIS Broach Size 7 Right 4845-2-958 ABGII MIS Broach Size 8 Right 4845-2-966 ABGII MIS Broach Size 6 Left 4845-2-967 ABGII MIS Broach Size 7 Left 4845-2-968 ABGII MIS Broach Size 8 Left 4845-2-951

ANATO Broach Size 1 Right

4845-2-952

ANATO Broach Size 2 Right

4845-2-953

ANATO Broach Size 3 Right

4845-2-954

ANATO Broach Size 4 Right

4845-2-955

ANATO Broach Size 5 Right

4845-7-956


ANATO Broach Size 6 Right

4845-7-957

ANATO Broach Size 7 Right

4845-7-958

ANATO Broach Size 8 Right



1020-1460 — Offset Broach Handle

#### Optional

1440-1460 Straight Broach Handle

NOTES

NOTES



325 Corporate Drive Mahwah, NJ 07430 t: 201 831 5000

stryker.com

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker's product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your sales representative if you have questions about the availability of products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: ABG, Anato, BixCut, Howmedica, Osteonics, Stryker, Unitrax, V40. All other trademarks are trademarks of their respective owners or holders.

BIOLOX delta is a registered trademark of Ceramtec Ag.

ANATO-SP-2\_19542

Copyright © 2019 Stryker Printed in USA