


# Hoffmann® 3

## External Fixation System




**Ready for all.**  
Yielding for nothing.

## Hoffmann 3 External Fixation System



### Delta couplings

- "Delta" design allows for 6 different pin/bar diameter connections through a single coupling
  - **Green:** Ø5mm, Ø8mm, Ø11mm
  - **Grey:** Ø4mm, Ø5mm, Ø6mm
- Inverted pin-to-rod offering for accessibility to thumbwheel in different frame configurations.
- Built-in, removable thumbwheels for provisional tightening
- Multiplanar Delta Couplings allow for increased flexibility
- Snap-fit technology designed for a fast ex-fix application
- Prevent "pop-out" with two turns of the thumbwheel



### Vectran coated carbon fiber rods

- Ø11mm carbon fiber rods
- 100mm – 650mm in length
- Vectran coating for electric insulation

### Pin clamps and rod couplers

- 5-hole and 10-hole pin clamps with customizable, fixed posts
- Modular or pre-welded versions available
- Straight, 30°, and 90° posts
- 30° fixed angle connection for Ø11mm rods
- Pelvic and knee bridging applications
- Asymmetric capability may aid in anatomic positioning



### Apex pin portfolio

- Self-drilling / self-tapping, blunt, cancellous, and transfixing Apex pins
- Stainless steel & titanium options, HA coating available
- Hybrid pins allow use of 3 & 4mm threads with the 5mm side of delta couplings
- Universal pin chuck allows stab-&-grab of all Apex pins



### Sterile field kits

- Field kits with pre-sterilized components and instruments needed for multiple indications
- Available for civilian use
- Bridging and emergency options available

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker's product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for Cleaning and Sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Apex, Hoffmann, Stryker. All other trademarks are trademarks of their respective owners or holders.

The products listed above are CE marked.

Content ID: H-SS-21, 05-2023  
Copyright © 2023 Stryker



#### Manufacturer:

Stryker GmbH  
Bohnackerweg 1  
2545 Selzach  
Switzerland

stryker.com