


Implant Extraction Set

System Components Guide
Module One & Two




System Components - Instruments



Catalog #	Description
1806-6150	Empty Tray Module 1
1806-6152	Module 1 (Complete)
Universal Screw Extraction Devices	
1806-6106	Screwdriver Bits Conical, Hex. 2.5mm
1806-6108	Screwdriver Bits Conical, Hex. 3.5mm
1806-6110	Screwdriver Bits Conical, Hex. 4.0mm
1806-6112	Screwdriver Bits, Hex. 5.0mm
1806-6114	Screwdriver Bits, Hex. 6.3mm
1806-6116	Screwdriver Bits, Hex. 8.0mm
1806-6118	Screwdriver Bits, Torx T 8.0
1806-6119	Screwdriver Bits, Torx T 15.0
1806-6120	Screwdriver Bits, Torx T 20.0
1806-6121	Screwdriver Bits, Torx T 25.0
1806-6113	Spreading Screwdriver Bit, Hex. 5.0mm
1806-6115	Spreading Screwdriver Bit, Hex. 6.3mm
1806-6117	Spreading Screwdriver Bit, Hex. 8.0mm
1806-6102	Teardrop Handle, Solid
1806-6104	Cannulated Teardrop Handle
1806-6104	Rotation Rod used with Cannulated Teardrop Handle
Universal Nail Extraction Device	
1806-0353	Conical Extraction Rod, Standard 6mm
1806-0350	Conical Extraction Rod, Standard 7/8mm
1806-6125	Conical Extraction Rod, 10mm
1806-6130	Conical Extraction Rod, Male, 13mm
1806-0130	Wrench 8/10mm, Standard
1806-0150	Strike Plate, Standard
1806-6100	Sliding Hammer
1806-0110	Universal Rod
1806-0170	Slotted Hammer, Standard
Universal Lag Screw Extraction Device	
1806-6140	Connector Gamma
1806-6141	Connector Gamma U-Blade
1806-6142	Connector Gamma3 U-Blade
1806-6144	Connector Gamma3
1806-6143	Connector Dyax (Asiatic)
1806-6146	Connector Omega
1806-6148	Connector OHS/OCS
1806-6136	Threaded Rod M7
1806-6137	Threaded Rod M5
1806-6138	Threaded Rod M4
1114-5002	Spanner 17mm Standard
1806-6135	Extraction Instrument Lag Screw
1806-6139	Nut

System Components - Instruments



Catalog # Description

Universal Emergency Screw Extraction Devices

1806-6151	Empty Tray Module 2
1806-6153	Module 2 (Complete)
1806-6169	Conical Extractor, Male, Small
1806-6170	Conical Extractor, Male, 2.5mm
1806-6171	Conical Extractor, Male, 3.5mm
1806-6172	Conical Extractor, Male, 4.0mm
1806-6173	Conical Extractor, Male, 5.0mm
1806-6165	Crown Drill (Trephine), 3.0mm
1806-6166	Crown Drill (Trephine), 4.0mm
1806-6167	Crown Drill (Trephine), 5.0mm
1806-6168	Crown Drill (Trephine), 6.5mm
1806-6180	Crown Drill (Trephine), 8.0mm
1806-6181	Crown Drill (Trephine), 10.0mm
1806-6182	Crown Drill (Trephine), 12.0mm
1806-6183	Conical Extractor, Female, 3.0mm
1806-6184	Conical Extractor, Female, 4.0mm
1806-6185	Conical Extractor, Female, 5.0mm
1806-6186	Conical Extractor, Female, 6.3mm
1806-6175	Drift Punch, 2.7mm
1806-6176	Drift Punch, 3.7mm
1806-6177	Drift Punch, 5.0mm
1806-6102	Teardrop Handle, Solid
1806-6160	Extraction Hook (can also be stored in Module 1), Small
1806-6162	Extraction Hook (can also be stored in Module 1), Large
1806-6178	Forceps for Screw Removal, Small
1806-6179	Forceps for Screw Removal, Large

Reconstructive

Hips
Knees
Trauma & Extremities
Foot & Ankle
Joint Preservation
Orthobiologics & Biosurgery

MedSurg

Power Tools & Surgical Accessories
Computer Assisted Surgery
Endoscopic Surgical Solutions
Integrated Communications
Beds, Stretchers & EMS
Reprocessing & Remanufacturing

Neurotechnology & Spine

Craniomaxillofacial
Interventional Spine
Neurosurgical, Spine & ENT
Neurovascular
Spinal Implants



For detailed information concerning the indications of the product(s) please refer to the dedicated Operative Technique.

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate a Stryker product.

A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for Cleaning and Sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders. The above listed products are CE Marked.

IES-BR-1 Rev 1, 05-2014

Copyright © 2014 Stryker

Manufactured by:

Stryker Trauma GmbH
Prof.-Küntschers-Straße 1–5
D - 24232 Schönkirchen
Germany

www.stryker.com