

# Hoffmann® LRF

## Circular External Fixation System

**Stable.**


Radiolucent carbon fiber rings offer enhanced stiffness compared to market-leading aluminium and carbon fiber rings\*

**Smart.**

Built-in dynamization mechanisms offer the option of 0-5mm of controlled dynamization

**Simple.**

Unique, 1 piece pin clamp design allows for multiplanar snap-fit connection of 3, 4, 5 and 6mm diameter Apex Pins



Stryker quality. Stryker experience. **Stryker Trauma & Extremities.**

Learn more about these and other Stryker products today. Visit [www.stryker.com](http://www.stryker.com).

\* Comprehensive Report of the Verification Tests Performed on the Hoffmann LRF System 2013, pages 22-25 (BML 13-049).

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Apex, Hoffmann, Stryker. All other trademarks are trademarks of their respective owners or holders.

H-FL-6\_Rev-2, 09-2020 • Copyright © 2016 Stryker.