
stryker

Allograft Wedges foot and ankle correction



stryker

Allograft wedges


foot and ankle correction

Cotton wedge

Ref #	Α	В	С
3102-1900	6mm	l4mm	24mm

Evans wedges

Ref #	Α	В	С
3102-1908	8mm		
3102-1910	10mm	20mm	22mm
3102-1912	12mm		



Wedge trials

Wedge trials allow surgeons to efficiently determine ideal wedge size intraoperatively. The trials were designed with a convenient full-wedge or height-only sizing option with heights of 6mm, 8mm, 10mm, and 12mm.

Ref #	Description	
705760	Wedge trial, 6mm	
705761	Wedge trial, 8mm	
705762	Wedge trial, 10mm	
705763	Wedge trial, 12mm	
940472	Wedge trial insert for VariAx 2 modular tray	





Trauma & Extremities

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for Cleaning and Sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, VariAx. All other trademarks are trademarks of their respective owners or holders.

Content ID: AGFT-BR-4 Rev-1, 04-2022 Copyright © 2022 Stryker



Wedge trials manufactured by:

Stryker 325 Corporate Dr. Mahwah, NJ 07430



Allograft wedges manufactured by:

AlloSource 6278 S Troy Cir. Centennial, CO 80111 t: 800 557 3587 stryker.com