

Trauma & Extremities

SOMA® designed



Implants Before SOMA

- Before SOMA, implants were designed using bones from cadavers, limited 2D X-ray images, and subjective design input.1
- Minimal 3-dimensional data limited the engineers' ability to extrapolate their data to the whole population.2
- Achieving true morphological fit in Orthopaedic prostheses has been an elusive industry goal.3

Designing with SOMA

- SOMA consists of a database of over 19,500 3D bones generated from CT scans of patients. This comprehensive database and proprietary software helps design a better fitting implant and allows for true evidence based design.1
- The 3D computer-aided design modeling continuously provides virtual feedback and prevents testing of nonconforming prototypes thereby enabling a more efficient product development process.2
- SOMA was first used to design T2 GTN in 2008, and since then has contributed to over 20 innovative Stryker products.

Innovative and differentiated

Understanding the technology

- The proprietary software to evaluate bone morphology and anatomical fit consists of three components
 - The Bone database management tool (BODAMAT) selects a target population based on specifics entered (i.e. age, gender, ethnicity, etc.)2
 - Stryker Anatomy Analysis Tool (SAAT) allows for quantitative analyzation of shape variation of bones and aids in identifying geometric measurements critical for design.2
 - **Stryker Implant Fitting Tool (SIFT)** performs an analysis of how well an implant fits on the bones in the SOMA database and enables comparison of our plate fit to competitive plates.2

Why fit matters

- A well-fitting plate may reduce the need for bending during the procedure which is important for preserving locking technology.2
- Reduced need for intra operative plate bending may allow for improved OR efficiency.2
- A well-fitting plate may result in reduction of soft tissue impingement and might decrease the risk of skin irritations.2
- A well-fitting plate can help with fracture alignment.2

stryker

PRO System



VariAx Clavicle



AxSOS 3 Titanium

Locking Plate System



Reasons to believe

- SOMA's proprietary software enables evidenced based design and the ability to verify fit.^{2,4}
- Recent plates designed using SOMA include VariAx Clavicle, PRO pelvis plates, and AxSOS 3 distal anterolateral tibia, proximal medial tibia, and the distal medial tibia plates.
- Enhanced anatomical fit in SOMA designed plates has been validated by feedback from our customers.⁵
- Stryker continues to innovate using SOMA. Several new tools, such as bone density analysis, are currently being developed to aid in future implant designs.³

"Using an underlying database with accompanying computation tools such as **SOMA** can be a powerful and efficient approach towards the development and advancement of osteosynthesis plates in trauma surgery, ultimately resulting in plates with high levels of **anatomic compliance** and **potential clinical benefits.**"²

References

- 1. Medical Education External SOMA Presentation, SOMA-EM-9_18207
- Petersik, A., Homeier, A., Hoare, S. G., Oldenburg, G. V., Gottschling, H., Schröder, M., & Burgkart, R. (2018). A numeric approach for anatomic plate design. Injury, 49, S96-S101.
- 3. Stryker Orthopaedics and Analytics (SOMA): A Review, SOMA-AR-2,01-2018
- 4. Competitive Fitting Study for Locking Proximal lateral Humerus Plates, AxSOS-WP-6 Rev 1
- 5. SOMA-EM-4, 12-2016

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: AxSOS, PRO, SOMA, Stryker, Stryker Orthopaedics, VariAx. All other trademarks are trademarks of their respective owners or holders.

Manufactured by:

Stryker GmbH Bohnackerweg 1 2545 Selzach Switzerland

stryker.com