


T2 ICF

Intramedullary nailing system

Nailing down complex cases.



T2 ICF

Intramedullary nailing system

stryker

Strong


- ▶ Designed to allow up to 10mm of internal compression⁵
- ▶ Large nail diameters (8, 9, 10mm) designed to help prevent nail breakage^{1,2}
- ▶ Type II Anodized titanium nail and endcap¹



Diameter	Length (10mm increments)
Ø8.0mm	80 to 140mm
Ø9.0mm	80 to 140mm
Ø10.0mm	80 to 140mm

Stable

- ▶ Intended for foot stabilization and alignment
- ▶ Nail designed to be compatible with cross screws to help prevent nail back-out²
- ▶ Nail and endcap design via SOMA⁶



*4x higher fatigue strength compared to predicate [501±14 N vs. 121±14 N, p=0.001]¹
 **2x higher static axial load to failure than predicate [2859±180N vs. 1204±258N, p<0.001]²


Efficient

- ▶ System instrumentation designed to promote:
 - An efficient and reproducible procedure^{3,7}
 - Accurate nail insertion^{3,4}
- ▶ Radiolucent instrumentation for clear visualization of implant positioning^{3,5}

New instruments



Guide Wire Targeter



Targeting Arm

References:

1. Stryker V18359 (A0045318, Ver. 1 7/1/2019) T2 ICF Implant / Synthes Midfoot Fusion Bolt Cantilever Bending Test; 2. Stryker V18361 (A0045344, Ver. 2 3/2/2020) T2 ICF Static axial load to failure; 3. Stryker V18368 (A0050277, Ver. 1 2/25/2020) T2 ICF Targeting System – System Handling Test; 4. Stryker V18371 (A0046580, Ver. 3 2/28/2020) T2 ICF Targeting System_ Life Cycle Simulation & Stiffness Test; 5. Stryker T2 ICF Op Tech T2-ST-32,07-2019; 6. Stryker A0043891 - T2ICF-DOF 30-502 SOMA analysis report- SAAT_2019_04 1.1; 7. Stryker V19379 Ver. 2 (A0051245) T2ICF Guide Wire Targeter-Shuttle-Final Design

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, T2. All other trademarks are trademarks of their respective owners or holders.

This document is only approved for use in the USA.

T2-SS-10, 06-2020

Copyright © 2020 Stryker

Manufacturer:

Stryker GmbH
 Bohnackerweg 1
 2545 Selzach,
 Switzerland

stryker.com