


**Spine** 

## AlloCraft® CL Cervical Allograft Spacer



- CorticoCancellous
- 5° of lordosis
- Chamfered sides
- Freeze dried

## AlloCraft CL

AlloCraft CL is a cervical allograft spacer designed for Anterior Cervical Discectomy and Fusion (ACDF) procedures. This cervical spacer features a cortical block to provide stability and a cancellous center, which may help serve as a scaffold for new bone growth. The graft is freeze dried and has four chamfered edges to ease insertion.



Part Number	Description
6183-5-005	5mm AlloCraft CL Spacer
6183-5-006	6mm AlloCraft CL Spacer
6183-5-007	7mm AlloCraft CL Spacer
6183-5-008	8mm AlloCraft CL Spacer
6183-5-009	9mm AlloCraft CL Spacer
6183-5-010	10mm AlloCraft CL Spacer

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: AlloCraft, Stryker. All other trademarks are trademarks of their respective owners or holders.

BIALL-SS-1 SC/GS 04/14 2 Pearl Court Allendale, NJ 07401 t: 201-760-8000

www.stryker.com