
<u>stryker</u>

Spine

AVS Anchor®-C Cervical Cage System


- Large Graft Volume
- Angled and Flexible Instrumentation
- Self-Locking Screws



AVS Anchor®-C Cervical Cage System

Optimized Screw Angulation

35° cephalad/caudal 10° medial convergence



Self-Locking Screws


Fewer surgical steps may allow for a faster procedure



Fixed Angle Screws

Self-tapping and Self-drilling Ø3.5mm: 8-14mm Ø4.0mm: 8-14mm





No Anterior Profile

Large Graft Volume

Designed for maximum bone graft-to-host bone contact



Adaptive Sizing

Heights: 6-12mm 0°, 4°, and 8° Lordosis:

Footprints: 12x14mm and 14x16mm

Intuitive Instrumentation

Angled, flexible, and all-in-one instrument options designed for:

- Potentially faster surgery due to fewer instruments needed
- Easier access to more challenging levels of the cervical spine
- Visual screw locking confirmation



A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: AVS, Stryker. All other trademarks are trademarks of their respective owners or holders.

CVANC-SS-1

SC/GS 07/14

EU Operations

ZI Marticot 33610 Cestas - France t: +33 (0) 5 57 97 06 30 f: +33 (0) 5 57 97 06 31

Allendale, New Jersey 07401

+1 201 760 8108

www.stryker.com

Phone: +1 201 760 8000

US Operations

2 Pearl Court

Fax:

Web: