

AlloWrap® DS Amniotic Membrane



- Versatile
- Dual layer
- **Ease of use**



AlloWrap DS

Amniotic Membrane


Amniotic membranes are harvested from a growth factor rich environment. They are commonly used in a variety of applications, including wound coverage, as a barrier and for nerve or tendon protection.¹

Dual layer

- The double-sided membrane features 2 layers of amniotic tissue, oriented with the epithelial layers facing outward.
- This processing allows for omni-directional implantation



- No product preparation is required
- Product is 0.1–0.2mm thick
- Room temperature storage
- Terminally sterilized 10⁻⁶



Ref#	Product description	
3102-2002	AlloWrap DS 2 x 2cm	
3102-2004	AlloWrap DS 2 x 4cm	
3102-2006	AlloWrap DS 4 x 4cm	
3102-2008	AlloWrap DS 4 x 8cm	

AlloWrap DS 2 x 2cm 3102-2002	AlloWrap DS - 2 x 4cm 3102-2004	AlloWrap DS	AlloWrap DS
AlloWrap DS 2 x 2cm 3102-2002		4 x 4cm 3102-2006	4 x 8cm 3102-2006

References:

 Duarte IGL et al. Amniotic Membrane as a biological dressing in infected wound healing in rabbits. Acta Cirurgica Brasileria (2014). Vol. 29 (4) p334-339.

Trauma & Extremities

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for Cleaning and Sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders. AlloWrap is a registered trademark of and is manufactured by AlloSource.



Manufactured by:

AlloSource 6278 S. Troy Cir. Centennial, CO 80111 t: 800 557 3587

stryker.com