




# AVS Anchor®-L

## Lumbar Cage System

- One Step Locking Plate
- Variable Screw Angulation
- Simple Instrumentation




**Integrated Fixation**



**Supplemental Fixation**


# AVS Anchor-L Lumbar Cage System

## One Step Locking Plate



- Reduces locking steps to one
- Absence of screw lock enables lagging for increased compression of the graft
- Simple Plate Inserter

## Variable Screw Angulation



- Allows 35 - 45° of sagittal angulation
- Self-drilling, self-tapping, dual lead screws
- Tapered minor diameter for increased insertion torque<sup>1</sup>

## Simple Instrumentation



- Standard, Flexible Shaft, U-Joint, Ball Hex and 90° Screwdriver options
- Simple All-in-One Guide with removable shaft
- Same hex for entire system

## One Device – Two Ways to Stabilize

The latest in innovation from **Stryker Spine**

<sup>1</sup> Inceogiu,S., Ferrara, L. McClain, R. (2004) Pedicle Screw Fixation Strength: pullout versus insertional torque. The Spine Journal, 4, 513-518.

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: AVS Anchor, Stryker. All other trademarks are trademarks of their respective owners or holders.

Literature Number: TLANLSS13060\_GLB  
SC/GS 07/13

Copyright © 2013 Stryker  
Printed in USA

**US Operations**  
**Stryker Spine**  
2 Pearl Court  
Allendale, New Jersey 07401  
**Phone:** +1 201 760 8000  
**Fax:** +1 201 760 8108  
**Web:** [www.stryker.com](http://www.stryker.com)

**EU Operations**  
**Stryker Spine SAS**  
ZI Marticot  
33610 Cestas – France  
**t:** +33 (0) 5 57 97 06 30  
**f:** +33 (0) 5 57 97 06 31