

# Redesigning a K-Wire limits inadvertent advancement: A biomechanical study

Carolyn Hardin • Joshua E Preiss BS • Sergio Gutierrez PhD • Robert E Isaacs MD

## Introduction

The migration of K-Wires is a well-known complication that can occur with potentially catastrophic consequences. This study investigates a guidewire technology designed to increase the force it takes to advance a K-Wire after it is deployed.

## Methods

A biomechanical study was designed to test the force required to displace the guidewire through both cancellous and cortical bone. Guidewires were placed at every vertebral level between Ll to L5 in human cadaveric specimens. The design was internally randomized between a standard guidewire and a Y-Wire, allowing each vertebra to serve as its own control. Displacement and axial load were measured. The regions of interest were the mid vertebral cancellous bone and the anterior cortical wall.




Figure 1. Radiographic image of splayed Y-Wire in mid-vertebral cancellous bone

All authors: Duke University Medical Center Department of Surgery Division of Neurosurgery

#### Results

The force for a standard wire stayed roughly the same throughout the cancellous bone tested (p=0.83). In contrast, the force to advance a Y-Wire increased in a linear fashion. The median cancellous force values exceeded 400% of the values obtained from a standard guidewire (435%, range 272% - 2136%). The force to displace a Y-Wire exceeded the cortical values obtained in the same vertebra after it was deployed 5.0mm (range 1.0-8.4mm). Additionally, the force required to breach the anterior cortex with a Y-Wire was 219% greater than a standard guidewire (p<0.001, range 184%-7129%).




Figure 2. Average force for a guidewire to break through cortical Bone

#### Conclusions

Y-Wire offers the ability to increase the push-in force to advance a guidewire placed into cancellous bone to values that exceed the force which cause standard guidewires to break through the anterior cortex. By altering the distal tip design of a guidewire, Y-Wire substantially decreases the chances of inadvertently advancing a guidewire away from its desired location.

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

MITLF-BM-1 SC/GS 03/16

Copyright © 2016 Stryker Printed in USA



Stryker's Spine Division 2 Pearl Court Allendale, NJ 07401-1677 USA t: 201 749 8000 www.stryker.com